## HEPES缓冲液（（BZ228）500ml

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[HEPES buffer\(\(BZ228\)](#)

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Chemazone482

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4743

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[Biochemazone™](#)

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1000mL, 100mL, 200mL, 500mL

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Product name HEPES buffer (1×with Calcium-Magnesium) HEPES buffer (1×, containing calcium and magnesium)

Product specifications 500mL

Reagent storage Transport at room temperature, store at 2-8°C

Validity period: 12 months (use as soon as possible after opening the lid)

### Product Introduction

HEPES is 4-(2-Hydroxyethyl)piperazine-1-ethanesulfonic acid, which is non-toxic to cells. HEPES is a non-ionic amphoteric buffer, which can effectively control the pH in the range of 6.8 to 8.2, especially at pH 7.2 to 7.4. It has good buffering capacity. The final concentration is generally 10 to 50 mmol/L. mmol/L HEPES can achieve better buffer capacity.

Phygene HEPES buffer (1×, containing calcium and magnesium) is a commonly used solution in cell-cell adhesion experiments. It is mainly composed of HEPES, sodium chloride, potassium chloride, phosphate, glucose, etc., and contains Ca<sup>2+</sup>, Mg<sup>2+</sup>, pH7.4, sterilized by filtration. It is often used to clean tissues or cells (especially embryonic cells) during long-term aggregation culture to protect the cadherin of cells without affecting the formation of cell aggregation.

This product is sterilized by 0.2µm filter

### Use reference

1. Mainly used for washing cells, similar to PBS for direct use.
2. If it is a cell, just wash it directly. If it is a tissue, the tissue should be cut out, placed in pre-cooled HEPES buffer (1 ×, containing calcium and magnesium), and cut into small pieces, collected with the minimum amount of buffer, and transferred to a BSA-coated petri dish. In, use scissors to cut the tissues for downstream cell culture experiments.

### Precautions

1. Pay attention to aseptic operation and try to avoid contamination.
2. For your safety and health, please wear lab coats and disposable gloves.

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**HEPES buffer(1×with Calcium-Magnesium) and customer's specified configurations can be delivered upon request. Stabilizer, pH Value, Pack Size, and Contents of HEPES**

buffer(1xwith Calcium-Magnesium) are customizable.

Contact us for Customization: [sales@coreab.cn](mailto:sales@coreab.cn)

Transportation and storage: Transportation at room temperature, storage at 4 °C. For Longterm Storage Freeze it at -20 °C

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