## **PBST缓冲液 (1x, PH7.4) (BZ218) 1000ml**

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[PBST \(1x,pH7.4\) \(BZ218\)](#)

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Chemazone507

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8868

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Store at 2-8°C

## Product Introduction

Phosphate buffered saline, Phosphate-Buffered Saline, is one of the most commonly used phosphate buffered solutions, also known as Dulbecco's phosphate buffered saline, Dulbecco's PBS, D-PBSA, CMF-DPBS, D-PBS, DPBS.

PBST (1x, pH7.4), also known as PBS-Tween, is composed of sodium chloride, disodium hydrogen phosphate, sodium dihydrogen phosphate, etc., contains 0.05% Tween20, does not contain Ca<sup>2+</sup>, Mg<sup>2+</sup>, and has a pH of 7.4. It is mostly used ELISA cleaning solution can also be used to clean PVDF membrane or NC membrane in Western Blot experiment.

This product is 1x working fluid and can be used directly. Storage at 4°C can extend the shelf life of the product. Low temperature storage is not recommended for long-term storage.

## Instructions for use

1. Operation refer to the specific requirements of the experiment. No need to prepare, use directly.
2. When cleaning the membrane, repeat 3 times for 5-10 minutes each time.

## Precautions

1. Low temperature storage is easier to produce precipitation. If there is precipitation, please shake or place it in a 37°C water bath until it is completely dissolved before use.
2. For your safety and health, please wear lab coats and disposable gloves for operation.

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**PBST (1x,pH7.4) Buffer and customer's specified configurations can be delivered upon request. Stabilizer, pH Value, Pack Size, and Contents of PBST (1x,pH7.4)Buffer are customizable.**

**Contact us for Customization: [sales@coreab.cn](mailto:sales@coreab.cn)**

**Transportation and storage: Transportation at room temperature, storage at 4 °C. For Longterm Storage Freeze it at -20 °C**

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