CE Marked Material
Monitor Sample for Anti HTLV-I
NIBSC code: 03/104-xxx
Instructions for use
(Version 18.0, Dated 13/07/2023)

This material is an 'Annex II List A' IVD and complies with the requirements of the "EU in vitro diagnostic medical device directive 98/79/EC"

1. INTENDED USE

This product is CE marked for use as an IVD within the UK, EU member states and EEA countries. In all other territories this product can be used for research purposes only.

This product is CE marked for use as an IVD within the UK, EU member states and EEA countries. In all other territories this product can be used for research purposes only.

The Monitor Sample for Anti-HTLV-I manufactured by NIBSC has been on the market since 2003. It is intended for use in the field of in vitro diagnostics, in conjunction with diagnostic immunoassay test kits / systems for the detection of antibodies to HTLV-I to monitor the performance of these systems. It can be used to monitor the consistency of test performance using statistical process control on a daily basis and, over a period of time, as a retrospective monitor of batch performance. However it is for the user to establish suitability of purpose. It is expected that the Monitor Sample for Anti-HTLV-I will be detected by assays with high analytical sensitivity. The Monitor Sample for Anti-HTLV-I is for use in immunoassay test kits or systems that are capable of detecting antibodies to HTLV-I.

2. CAUTION

This preparation is not for administration to humans or animals in the human food chain.

This preparation contains material of human origin. It contains a dilution of a plasma pool prepared from 9 individual human plasma samples each known to be positive for antibodies to HTLV-I. The plasma pool and each plasma donation has been tested and found negative for HBsAg, anti-HCV and anti-HIV-1/2. Furthermore, the plasma pool has also been tested and found negative for HCV RNA. As with all materials of biological origin, this preparation should be regarded as potentially hazardous to health. It should be used and discarded according to your own laboratory's safety procedures. Such safety procedures should include the wearing of protective gloves and avoiding the generation of aerosols. Care should be exercised in opening ampoules or vials, to avoid cuts.

3. UNITAGE

No unitage is assigned to this material.

4. CONTENTS

Country of origin of biological material: United Kingdom.
Each vial of the Monitor Sample for Anti-HTLV-I contains approximately 3ml of liquid and consists of:-

- o Pooled and diluted anti-HTLV-I-positive human plasma.
- o Phosphate buffered saline (PBS-A) .
- o 5% Bovine serum albumin.
- o 0.05% Bronidox as preservative.

5. STORAGE

Upon receipt, vials should be stored at +2-8 C. Reference materials are held at NIBSC within assured, temperature-controlled storage facilities.

6. DIRECTIONS FOR OPENING

Vials have a screw cap; an internal stopper may also be present. The cap should be removed by turning anti-clockwise. Care should be taken to prevent loss of the contents. Please note: If a stopper is present on removal of the cap, the stopper should remain in the vial or be removed with the cap.

7. USE OF MATERIAL

The user should follow the assay kit manufacturer's instructions. The Monitor Sample should be used in addition to the positive and negative controls supplied with each kit and treated as if it were a routine test sample. Only the controls supplied with each kit by the manufacturer are used to determine the validity of assays and to calculate the cut-off value that forms the basis for donor screening or diagnosis.

THIS PREPARATION IS SUPPLIED FOR USE IN ITS FINAL FORM AND MUST NOT BE FURTHER DILUTED OTHER THAN AS REQUIRED IN INDIVIDUAL TEST PROCEDURES.

The Result Reporting System (RRS) has been developed by the National Institute for Biological Standards and Control (NIBSC) for the data monitoring of its serology and NAT quality control (QC) reagents. These include the Quality Control Reagent Unit (QCRU) and Clinical Virology Network (CVN) reagents. The system has been successfully running for serology assays for several years, collecting thousands of data points a year; and more recently for data derived from reagents used in Nucleic Acid-based Technologies (NAT) assays. <https://www.nibsc.org/products/rrs.aspx>

Users should be aware that different batches of the same assay kit and / or different assay conditions may give different values with this preparation.

8. STABILITY

The expiry date is stated on the vial label. Once opened, this preparation may be used for up to 3 months (as long as the expiry date is not exceeded), and should be stored at +2-8 C between each use.

The stability of this material is monitored by NIBSC. The stability of previous batches of the Monitor Sample for Anti-HTLV-I have not lost potency when stored at +2-8 C for periods of up to 36 months. The shelf life assigned to this preparation is 36 months when stored at +2-8 C.

Users who have data supporting any deterioration in the characteristics of any Working Standard / Monitor Sample are encouraged to contact NIBSC.

9. REFERENCES

Not applicable.

10. ACKNOWLEDGEMENTS

EC REP: Advena Ltd. Tower Business Centre, 2nd Flr., Tower Street, Swatar, BKR 4013 Malta

11. FURTHER INFORMATION

Further information can be obtained as follows;
 This material: enquiries@nibsc.org
 WHO Biological Standards:
<http://www.who.int/biologicals/en/>
 JCTLM Higher order reference materials:
<http://www.bipm.org/en/committees/jc/jctlm/>
 Derivation of International Units:
http://www.nibsc.org/standardisation/international_standards.aspx
 Ordering standards from NIBSC:
<http://www.nibsc.org/products/ordering.aspx>
 NIBSC Terms & Conditions:
http://www.nibsc.org/terms_and_conditions.aspx

12. CUSTOMER FEEDBACK

Customers are encouraged to provide feedback on the suitability or use of the material provided or other aspects of our service. Please send any comments to enquiries@nibsc.org

13. CITATION

In all publications, including data sheets, in which this material is referenced, it is important that the preparation's title, its status, the NIBSC code number, and the name and address of NIBSC are cited and cited correctly.

14. MATERIAL SAFETY SHEET

Classification in accordance with Directive 2000/54/EC, Regulation (EC) No 1272/2008: Not applicable or not classified

Physical and Chemical properties	
Physical appearance: Liquid	Corrosive: No
Stable: Yes	Oxidising: No
Hygroscopic: No	Irritant: Yes (contains Bronidox)
Flammable: No	Handling: See caution, Section 2
Other (specify):	Contains material of human origin.
Toxicological properties	
Effects of inhalation: Not established. Avoid inhalation. May contain infectious HTLV-I	
Effects of ingestion: Not established. Avoid ingestion. May contain infectious HTLV-I	
Effects of skin absorption: Not established. Avoid contact with skin. May contain infectious HTLV-I	
Suggested First Aid	
Inhalation: Seek medical advice.	
Ingestion: Seek medical advice.	
Contact with eyes: Wash with copious amounts of water. Seek medical advice.	
Contact with skin: Wash thoroughly with water.	

Action on Spillage and Method of Disposal
Spillage of vial contents should be taken up with absorbent material wetted with an appropriate disinfectant. Rinse area with an appropriate disinfectant followed by water. Absorbent materials used to treat spillage should be treated as biological waste.

15. LIABILITY AND LOSS

In the event that this document is translated into another language, the English language version shall prevail in the event of any inconsistencies between the documents.

Unless expressly stated otherwise by NIBSC, NIBSC's Standard Terms and Conditions for the Supply of Materials (available at www.nibsc.org/About_Us/Terms_and_Conditions.aspx or upon request by the Recipient) ("Conditions") apply to the exclusion of all other terms and are hereby incorporated into this document by reference. The Recipient's attention is drawn in particular to the provisions of clause 11 of the Conditions.

16. INFORMATION FOR CUSTOMS USE ONLY

Country of origin for customs purposes*: United Kingdom * Defined as the country where the goods have been produced and/or sufficiently processed to be classed as originating from the country of supply, for example a change of state such as freeze-drying.
Net weight: 3g per vial.
Toxicity Statement: Non-toxic.
Veterinary certificate or other statement if applicable.
Attached: No

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