DATASHEET For Research Use Only

REAGENT SARS-CoV-2 Nucleocapsid Protein (His-tagged)

REAGENT REFERENCE 101042

LOT NUMBER: lot 2

PROVIDED: 100µg at 1.1 mg/mL of purified protein HEPES buffered saline, 5% glycerol,

1 mM EDTA, pH 8.

DESCRIPTION: N-terminal mini-histidine-tagged fusion protein derived from SARS-CoV-2

nucleocapsid protein (SARS-CoV-2). NCBI Reference Sequence: Accession No. YP_009724397. Residues 9 – 426 of the sequence below (N-terminal histidine tag underlined) match amino acids 2- 419 of YP 009724397.

PROTEIN SEQUENCE;

MHHHHHHGSDNGPQNQRNAPRITFGGPSDSTGSNQNGERSGARSKQRRPQGLPNNTASWF TALTQHGKEDLKFPRGQGVPINTNSSPDDQIGYYRRATRRIRGGDGKMKDLSPRWYFYYLGT GPEAGLPYGANKDGIIWVATEGALNTPKDHIGTRNPANNAAIVLQLPQGTTLPKGFYAEGSRGG SQASSRSSSRSNSSRNSTPGSSRGTSPARMAGNGGDAALALLLLDRLNQLESKMSGKGQQ QQGQTVTKKSAAEASKKPRQKRTATKAYNVTQAFGRRGPEQTQGNFGDQELIRQGTDYKHW PQIAQFAPSASAFFGMSRIGMEVTPSGTWLTYTGAIKLDDKDPNFKDQVILLNKHIDAYKTFPT EPKKDKKKKADETQALPQRQKKQQTVTLLPAADLDDFSKQLQQSMSSADSTQA

PURITY: Estimated by SDS-PAGE Nucleocapsid with less than 5% E. coli

contaminant proteins.

MOLECULAR WEIGHT: Experimental by Electrospray Ionisation MS, M+H+ ion 46508.08 is within


0.003% of calculated mass of M+H+ ion (46506.60).

HOMOGENEITY: ~ 95% by SDS-PAGE. Observable bands identified by MS as the expected

protein.


ANALYTICAL DATA:

User Ref: 101042	Version: 1.0	Version Date: 08/02/2022
Latest version of document available at: Q-DOCS		Date Printed: 10/02/2022 11:55
Document ID: 9234	Page 1 of 3	Issue Status: Published



SDS-PAGE Analysis of His-tagged NCP

Protein loaded at 10, 3 or 1 μ g/ well analyzed on a 10% polyacrylamide gel, stained with Coomassie Blue. Proteomic analysis indicates *E. coli* protein contamination of less than 1% in main band **A**. Bands **B** and **C** identified as fragments of NCP by proteomic analysis.



Absorbance spectrum determined in aqueous buffer over wavelength range ~240 -320 nm

STORAGE: Keep at -80°C. Avoid freeze-thaw cycles as reagent degradation may result.

APPLICATION: Suitable for immunoassay.

DEPOSITOR: Prof Jon Sayers, The University of Sheffield, UK

ACKNOWLEDGMENTS

Acknowledgment for publications should read "The following reagent was obtained from the Centre For AIDS Reagents, NIBSC, UK: SARS-CoV-2 Nucleoprotein (His-tagged) (#101042) from Prof Jon Sayers".

User Ref: 101042	Version: 1.0	Version Date: 08/02/2022
Latest version of document available at: Q-DOCS		Date Printed: 10/02/2022 11:55
Document ID: 9234	Page 2 of 3	Issue Status: Published

MATERIAL SAFETY SHEET

Physical properties (at room temperature)				
Clear, liquid				
None				
Chemical properties				
Yes	Corrosive:	No		
No	Oxidising:	No		
No	Irritant:	No		
	Clear, liquid None Chem Yes No	Clear, liquid None Chemical properties Yes Corrosive: No Oxidising:		

Other:

This product is a recombinant protein; It is the responsibility of the end user to seek local biosafety approval for the storage and handling of the material in their workplace.

Handling:

CAUTION - This preparation is not for administration to humans or animals in the human food chain. As with all materials of biological origin, this preparation should be regarded as potentially hazardous to health. It should be used and discarded according to your own laboratory's safety procedures. Such safety procedures should include the wearing of protective gloves and avoiding the generation of aerosols.

Toxicological properties				
Effects of inhalation:	Not established, avoid inhalation			
Effects of ingestion:	Not established, avoid ingestion			
Effects of skin absorption:	Not established, avoid contact with skin			
Suggested First Aid				
Inhalation	Seek medical advice			
Ingestion	Seek medical advice			
Contact with eyes	Wash with copious amounts of water. Seek medical advice.			
Contact with skin	Wash thoroughly with water.			
	Action on Spillage and Method of Disposal			

User Ref: 101042	Version: 1.0	Version Date: 08/02/2022
Latest version of document available at: Q-DO	Date Printed: 10/02/2022 11:55	
Document ID: 9234	Page 3 of 3	Issue Status: Published

Absorbent materials used to treat spillage should be treated as biologically hazardous waste.

Spillage of vial contents should be taken up with absorbent material wetted with disinfectant. Rinse

area with disinfectant followed by water.