



**International Ref. Preparation  
1st WHO International Reference Preparation for HIV-1 CRF's  
NIBSC code: 13/214  
Instructions for use  
(Version 1.0, Dated 02/04/2014)**

**1. INTENDED USE**

It has been known for some time that different subtypes of HIV-1 exist. There is the major group known as M consisting of subtypes A –J and more diverse groups of outliers such as group N, P and O. Initially nucleic acid-based tests had a narrow band of specificity targeting the B clade as this was most predominate in Europe and the US. Improvements have since been made in assay design in an attempt to detect a wider range of subtypes. However, it is recognised that recombination events have lead to a wider diversity of HIV-1 subtypes and in recent years viruses containing one or more subtype sequence within the genome have become evident, these viruses have become known as Circulating Recombinant Forms (CRF's).

It is known from collaborative studies conducted using an HIV-1 NAT panel containing many common subtypes that predominately circulate within in Europe, North America and Asia that some assays are still poor at detecting such subtypes, thus giving a low or negative result on samples that are known to be positive. In order to allow manufacturers of assays and laboratories running in house assays to validate the assays ability to detect CRF's, the WHO endorsed the development of the 1st HIV-1 NAT CRF panel

**2. CAUTION**

**This preparation is not for administration to humans or animals in the human food chain.**

The preparation contains material of human origin, and either the final product or the source materials, from which it is derived, have been tested and found negative for HBsAg, anti-HIV and HCV RNA.

As with all materials of biological origin, this preparation should be regarded as potentially hazardous to health. It should be used and discarded according to your own laboratory's safety procedures. Such safety procedures should include the wearing of protective gloves and avoiding the generation of aerosols. Care should be exercised in opening ampoules or vials, to avoid cuts.

**3. UNITAGE**

Due to the variation seen in different assays ability at detecting some of these panel members, no unitage has been assigned. However in order to gauge how each panel member behaved in different assays users are directed to the WHO ECBS report as stated in the reference section

**4. CONTENTS**

Country of origin of biological material: United Kingdom.

The final material was formulated and produced in the UK, diluent plasma was from a UK source however inactivated viruses contained within the panel originated from several different countries

**5. STORAGE**

On receipt the panel should be placed at -20 until use.

**Please note: because of the inherent stability of lyophilized material, NIBSC may ship these materials at ambient temperature.**

**6. DIRECTIONS FOR OPENING**

Vials have a screw cap; an internal stopper may also be present. The cap should be removed by turning anti-clockwise. Care should be

taken to prevent loss of the contents. Please note: If a stopper is present on removal of the cap, the stopper should remain in the vial or be removed with the cap.

**7. USE OF MATERIAL**

No attempt should be made to weigh out any portion of the freeze-dried material prior to reconstitution

Each vial should be reconstituted with 1ml of molecular grade water just prior to use. Vials should be left for a period of 20 minutes with occasional agitation, a visual check should be made to ensure all contents have fully reconstituted prior to extraction.

The material should be assayed neat or at serial dilutions to determine assay specificity and sensitivity.

**8. STABILITY**

Reference materials are held at NIBSC within assured, temperature-controlled storage facilities. Reference Materials should be stored on receipt as indicated on the label.

This panel should be used immediately following successful reconstitution. The reconstituted material should not be subjected to further storage or freeze thawing.

**9. REFERENCES**

WHO ECBS report WHO/BS/2013.2226

**10. ACKNOWLEDGEMENTS**

We would like to acknowledge the collaborative study group that analysed this panel as part of the establishment of this material.

**11. FURTHER INFORMATION**

Further information can be obtained as follows;

This material: [enquiries@nibsc.org](mailto:enquiries@nibsc.org)

WHO Biological Standards:

<http://www.who.int/biologicals/en/>

JCTLM Higher order reference materials:

<http://www.bipm.org/en/committees/jc/jctlm/>

Derivation of International Units:

[http://www.nibsc.org/standardisation/international\\_standards.aspx](http://www.nibsc.org/standardisation/international_standards.aspx)

Ordering standards from NIBSC:

<http://www.nibsc.org/products/ordering.aspx>

NIBSC Terms & Conditions:

[http://www.nibsc.org/terms\\_and\\_conditions.aspx](http://www.nibsc.org/terms_and_conditions.aspx)

**12. CUSTOMER FEEDBACK**

Customers are encouraged to provide feedback on the suitability or use of the material provided or other aspects of our service. Please send any comments to [enquiries@nibsc.org](mailto:enquiries@nibsc.org)

**13. CITATION**

In all publications, including data sheets, in which this material is referenced, it is important that the preparation's title, its status, the NIBSC code number, and the name and address of NIBSC are cited and cited correctly.

**14. MATERIAL SAFETY SHEET**

Classification in accordance with Directive 2000/54/EC, Regulation (EC) No 1272/2008: Not applicable or not classified

Physical and Chemical properties		
Physical appearance: Lyophilised	Corrosive:	No
Stable: Yes	Oxidising:	No



Hygroscopic:	No	Irritant:	No
Flammable:	No	Handling:	See caution, Section 2
Other (specify):			
<b>Toxicological properties</b>			
Effects of inhalation:	Not established, avoid inhalation		
Effects of ingestion:	Not established, avoid ingestion		
Effects of skin absorption:	Not established, avoid contact with skin		
<b>Suggested First Aid</b>			
Inhalation:	Seek medical advice		
Ingestion:	Seek medical advice		
Contact with eyes:	Wash with copious amounts of water. Seek medical advice		
Contact with skin:	Wash thoroughly with water.		
<b>Action on Spillage and Method of Disposal</b>			
Spillage of ampoule contents should be taken up with absorbent material wetted with an appropriate disinfectant. Rinse area with an appropriate disinfectant followed by water. Absorbent materials used to treat spillage should be treated as biological waste.			

#### 15. LIABILITY AND LOSS

In the event that this document is translated into another language, the English language version shall prevail in the event of any inconsistencies between the documents.

Unless expressly stated otherwise by NIBSC, NIBSC's Standard Terms and Conditions for the Supply of Materials (available at [http://www.nibsc.org/About\\_Us/Terms\\_and\\_Conditions.aspx](http://www.nibsc.org/About_Us/Terms_and_Conditions.aspx) or upon request by the Recipient) ("Conditions") apply to the exclusion of all other terms and are hereby incorporated into this document by reference. The Recipient's attention is drawn in particular to the provisions of clause 11 of the Conditions.

#### 16. INFORMATION FOR CUSTOMS USE ONLY

<b>Country of origin for customs purposes*:</b> United Kingdom * Defined as the country where the goods have been produced and/or sufficiently processed to be classed as originating from the country of supply, for example a change of state such as freeze-drying.
<b>Net weight:</b> 1g
<b>Toxicity Statement:</b> Toxicity not assessed
<b>Veterinary certificate or other statement</b> if applicable. <b>Attached:</b> No