



**WHO International Standard  
Prostate specific antigen (human) (total: PSA-ACT + free PSA)  
NIBSC code: 17/100  
Instructions for use  
(Version 2.0, Dated 29/01/2019)**

### 1. INTENDED USE

This material replaces the 1st International Standard (IS) for Prostate Specific Antigen (PSA) (90:10), coded 96/670, which has been widely used for the calibration of immunoassays for total PSA. Stocks of the 1st IS are exhausted and the WHO Expert Committee on Biological Standardization (ECBS) recognized (2014) the need for a replacement IS. The 1st IS for PSA contained seminal plasma-derived PSA, 90% bound to  $\alpha$ 1-antichymotrysin (PSA-ACT) and 10% in the free form and had an assigned content of 1  $\mu$ g total PSA per vial (1). The 2nd IS prepared in ampoules coded 17/100, contains seminal-fluid derived PSA-ACT and non-complexing (free) PSA. The 2nd WHO IS for PSA (PSA-ACT + free PSA), 17/100, was established at the 69th Meeting of WHO ECBS (2018). This material replaces the 1st IS for PSA (90:10), 96/670, which is discontinued.

### 2. CAUTION

**This preparation is not for administration to humans or animals in the human food chain.**

The preparation contains material of human origin, and either the final product or the source materials, from which it is derived, have been tested and found negative for HBsAg, anti-HIV and HCV RNA. As with all materials of biological origin, this preparation should be regarded as potentially hazardous to health. It should be used and discarded according to your own laboratory's safety procedures. Such safety procedures should include the wearing of protective gloves and avoiding the generation of aerosols. Care should be exercised in opening ampoules or vials, to avoid cuts.

### 3. UNITAGE

Each ampoule of the International Standard contains 0.50  $\mu$ g PSA

### 4. CONTENTS

Country of origin of biological material: USA.

Each ampoule contains the residue, after freeze-drying, of a 1.0 ml volume of a solution of 20 mM sodium phosphate buffer pH 7.4, 150 mM sodium chloride, 10 g/L bovine serum albumin (BSA) and 0.50  $\mu$ g/ml total prostate specific antigen comprised of PSA-ACT and free PSA.

### 5. STORAGE

Ampoules should be stored at -20°C.

**Please note: because of the inherent stability of lyophilized material, NIBSC may ship these materials at ambient temperature.**

### 6. DIRECTIONS FOR OPENING

DIN ampoules have an 'easy-open' coloured stress point, where the narrow ampoule stem joins the wider ampoule body. Various types of ampoule breaker are available commercially. To open the ampoule, tap the ampoule gently to collect material at the bottom (labelled) end and follow manufacturer's instructions provided with the ampoule breaker.

### 7. USE OF MATERIAL

**No attempt should be made to weigh out any portion of the freeze-dried material prior to reconstitution**

For all practical purposes each ampoule contains the same quantity of the substances listed above. Depending on the intended use, dissolve the total contents of the ampoule in a known volume of a suitable diluent. Users should make their own investigations into the type of diluent suitable for their use. If extensive dilutions are prepared, a carrier protein

should be added. The ampoules do not contain bacteriostat and solutions of the material should not be assumed to be sterile.

### 8. STABILITY

Reference materials are held at NIBSC within assured, temperature-controlled storage facilities. Reference Materials should be stored on receipt as indicated on the label. Analysis of accelerated thermal degradation samples of 17/100, measured by participants in a collaborative study (2) showed a predicted loss of 0.028 % of total PSA immunoreactivity per year when stored at -20°C. NIBSC follows the policy of WHO with respect to its reference materials. It is the policy of WHO not to assign an expiry date to their international reference materials. They remain valid with the assigned potency and status until withdrawn or amended. Reference materials are held at NIBSC within assured, temperature-controlled storage facilities. Reference materials should be stored on receipt as indicated on the label. In addition, once reconstituted, diluted or aliquoted, users should determine the stability of the material according to their own method of preparation, storage and use. Users who have data supporting any deterioration in the characteristics of any reference preparation are encouraged to contact NIBSC.

### 9. REFERENCES

- (1) Rafferty, B., Rigsby, P., Rose, M., Stamey, T. and Gaines Das, R. (2000) Reference Reagents for Prostate-specific Antigen (PSA): Establishment of the First International Standards for Free PSA and PSA (90:10). Clin Chem 46, 1310-1317
- (2) Ferguson, J., Atkinson, E., Rigsby, P. and Burns C. (2018) WHO International Collaborative Study of the Proposed 2<sup>nd</sup> WHO International Standard for Total PSA (PSA-ACT + free PSA). [http://www.who.int/biologicals/BS.2018.2340\\_PSA\\_Total.pdf?ua=1](http://www.who.int/biologicals/BS.2018.2340_PSA_Total.pdf?ua=1)

### 10. ACKNOWLEDGEMENTS

We gratefully acknowledge the important contributions of all the participants in the collaborative study and Dina Patel, UK NEQAS Immunology, Immunochemistry & Allergy, Sheffield, UK.

### 11. FURTHER INFORMATION

Further information can be obtained as follows;  
This material: [enquiries@nibsc.org](mailto:enquiries@nibsc.org)  
WHO Biological Standards:  
<http://www.who.int/biologicals/en/>  
JCTLM Higher order reference materials:  
<http://www.bipm.org/en/committees/jc/jctlm/>  
Derivation of International Units:  
[http://www.nibsc.org/standardisation/international\\_standards.aspx](http://www.nibsc.org/standardisation/international_standards.aspx)  
Ordering standards from NIBSC:  
<http://www.nibsc.org/products/ordering.aspx>  
NIBSC Terms & Conditions:  
[http://www.nibsc.org/terms\\_and\\_conditions.aspx](http://www.nibsc.org/terms_and_conditions.aspx)

### 12. CUSTOMER FEEDBACK

Customers are encouraged to provide feedback on the suitability or use of the material provided or other aspects of our service. Please send any comments to [enquiries@nibsc.org](mailto:enquiries@nibsc.org)

### 13. CITATION

In all publications, including data sheets, in which this material is referenced, it is important that the preparation's title, its status, the NIBSC code number, and the name and address of NIBSC are cited and cited correctly.

### 14. MATERIAL SAFETY SHEET

Classification in accordance with Directive 2000/54/EC, Regulation (EC) No 1272/2008: Not applicable or not classified

<b>Physical and Chemical properties</b>	
Physical appearance: Freeze dried powder	Corrosive: No
Stable: Yes	Oxidising: No
Hygroscopic: Yes	Irritant: No
Flammable: No	Handling: See caution, Section 2
Other (specify): Contains material of human origin	
<b>Toxicological properties</b>	
Effects of inhalation:	Not established, avoid inhalation
Effects of ingestion:	Not established, avoid ingestion
Effects of skin absorption:	Not established, avoid contact with skin
<b>Suggested First Aid</b>	
Inhalation:	Seek medical advice
Ingestion:	Seek medical advice
Contact with eyes:	Wash with copious amounts of water. Seek medical advice
Contact with skin:	Wash thoroughly with water.
<b>Action on Spillage and Method of Disposal</b>	
Spillage of ampoule contents should be taken up with absorbent material wetted with an appropriate disinfectant. Rinse area with an appropriate disinfectant followed by water. Absorbent materials used to treat spillage should be treated as biological waste.	

#### 15. LIABILITY AND LOSS

In the event that this document is translated into another language, the English language version shall prevail in the event of any inconsistencies between the documents.

Unless expressly stated otherwise by NIBSC, NIBSC's Standard Terms and Conditions for the Supply of Materials (available at [http://www.nibsc.org/About\\_Us/Terms\\_and\\_Conditions.aspx](http://www.nibsc.org/About_Us/Terms_and_Conditions.aspx) or upon request by the Recipient) ("Conditions") apply to the exclusion of all other terms and are hereby incorporated into this document by reference. The Recipient's attention is drawn in particular to the provisions of clause 11 of the Conditions.

#### 16. INFORMATION FOR CUSTOMS USE ONLY

<b>Country of origin for customs purposes*:</b> United Kingdom * Defined as the country where the goods have been produced and/or sufficiently processed to be classed as originating from the country of supply, for example a change of state such as freeze-drying.
<b>Net weight:</b> 0.02 g
<b>Toxicity Statement:</b> Non-toxic
<b>Veterinary certificate or other statement</b> if applicable. <b>Attached:</b> No

#### 17. CERTIFICATE OF ANALYSIS

NIBSC does not provide a Certificate of Analysis for WHO Biological Reference Materials because they are internationally recognised primary reference materials fully described in the instructions for use. The reference materials are established according to the WHO Recommendations for the preparation, characterization and establishment of international and other biological reference standards [http://www.who.int/bloodproducts/publications/TRS932Annex2\\_Inter\\_bi\\_olefstandardsrev2004.pdf](http://www.who.int/bloodproducts/publications/TRS932Annex2_Inter_bi_olefstandardsrev2004.pdf) (revised 2004). They are officially endorsed by the WHO Expert Committee on Biological Standardization (ECBS) based on the report of the international collaborative study which established their suitability for the intended use.