



**WHO International Standard**  
**1st WHO International Standard for SARS-CoV-2 antigen**  
**NIBSC code: 21/368**  
**Instructions for use**  
**(Version 1.0, Dated 01/11/2022)**

§

**1. INTENDED USE**

This preparation contains 0.01% formaldehyde-inactivated SARS-CoV-2 Omicron (B.1.1.529, sub-variant BA.1). The reference is intended for the standardization, and evaluation of performance and sensitivity of SARS-CoV-2 antigen detection as well as for the calibration of secondary reference materials.

**2. CAUTION**

**The material is not of human or bovine origin. This preparation is not for administration to humans or animals**

**3. UNITAGE**

This material has an assigned unitage of 5000 International Units of SARS-CoV-2 antigen per ampoule.

**4. CONTENTS**

Country of origin of biological material: United Kingdom.  
Each ampoule is a lyophilate of a preparation that contained 0.25 mL of clarified SARS-CoV-2 culture supernatant of Omicron (B.1.1.529, sub-variant BA.1), inactivated with 0.01% formaldehyde and diluted ~1/15 in Copan Universal Transport Medium (UTM-RT).

**5. STORAGE**

This preparation should be stored at -20°C or below on receipt. Please note because of the inherent stability of lyophilized material, NIBSC may ship these materials at ambient temperature.

**6. DIRECTIONS FOR OPENING**

DIN ampoules have an 'easy-open' coloured stress point, where the narrow ampoule stem joins the wider ampoule body. Various types of ampoule breaker are available commercially. To open the ampoule, tap the ampoule gently to collect material at the bottom (labelled) end and follow manufactures instructions provided with the ampoule breaker.

**7. USE OF MATERIAL**

**No attempt should be made to weigh out any portion of the freeze-dried material prior to reconstitution**

This material is supplied lyophilized and before use should be reconstituted in 0.25 mL of ultra-pure water. Reconstituted material should be used on the day of reconstitution.

**8. STABILITY**

Reference materials are held at NIBSC within assured, temperature-controlled storage facilities. Reference Materials should be stored on receipt as indicated on the label.

NIBSC follows the policy of WHO with respect to its reference materials.

**9. REFERENCES**

National Institute for Biological Standards and Control,  
Potters Bar, Hertfordshire, EN6 3QG. T +44 (0)1707 641000, [nibsc.org](http://nibsc.org)  
WHO International Laboratory for Biological Standards,  
UK Official Medicines Control Laboratory

Fryer JF, Hockley JG, Rigsby P, Morris CL, and the Collaborative Study Group. A Collaborative Study to Evaluate the Proposed 1st WHO International Standard for SARS-CoV-2 Antigen WHO/BS/2022.2426. World Health Organization, Geneva, Switzerland.

**10. ACKNOWLEDGEMENTS**

We gratefully acknowledge the important contributions of the collaborative study participants. Funding for the development of this standard was provided by the WHO Health Emergencies Programme.

**11. FURTHER INFORMATION**

Further information can be obtained as follows;  
This material: [enquiries@nibsc.org](mailto:enquiries@nibsc.org)  
WHO Biological Standards:  
<http://www.who.int/biologicals/en/>  
JCTLM Higher order reference materials:  
<http://www.bipm.org/en/committees/jc/jctlm/>  
Derivation of International Units:  
[http://www.nibsc.org/standardisation/international\\_standards.aspx](http://www.nibsc.org/standardisation/international_standards.aspx)  
Ordering standards from NIBSC:  
<http://www.nibsc.org/products/ordering.aspx>  
NIBSC Terms & Conditions:  
[http://www.nibsc.org/terms\\_and\\_conditions.aspx](http://www.nibsc.org/terms_and_conditions.aspx)

**12. CUSTOMER FEEDBACK**

Customers are encouraged to provide feedback on the suitability or use of the material provided or other aspects of our service. Please send any comments to [enquiries@nibsc.org](mailto:enquiries@nibsc.org)

**13. CITATION**

In all publications, including data sheets, in which this material is referenced, it is important that the preparation's title, its status, the NIBSC code number, and the name and address of NIBSC are cited and cited correctly.

**14. MATERIAL SAFETY SHEET**

Classification in accordance with Directive 2000/54/EC, Regulation (EC) No 1272/2008: Not applicable or not classified

Physical and Chemical properties	
Physical appearance: Lyophilized	Corrosive: No
Stable: Yes	Oxidising: No
Hygroscopic: No	Irritant: No
Flammable: No	Handling: See caution, Section 2
Other (specify): N/A	
Toxicological properties	
Effects of inhalation:	Not established, avoid inhalation
Effects of ingestion:	Not established, avoid ingestion
Effects of skin absorption:	Not established, avoid contact with skin



<b>Suggested First Aid</b>	
Inhalation:	Seek medical advice
Ingestion:	Seek medical advice
Contact with eyes:	Wash with copious amounts of water. Seek medical advice
Contact with skin:	Wash thoroughly with water.
<b>Action on Spillage and Method of Disposal</b>	
Spillage of ampoule contents should be taken up with absorbent material wetted with an appropriate disinfectant. Rinse area with an appropriate disinfectant followed by water. Absorbent materials used to treat spillage should be treated as biological waste.	

#### 15. LIABILITY AND LOSS

In the event that this document is translated into another language, the English language version shall prevail in the event of any inconsistencies between the documents.

Unless expressly stated otherwise by NIBSC, NIBSC's Standard Terms and Conditions for the Supply of Materials (available at [http://www.nibsc.org/About\\_Us/Terms\\_and\\_Conditions.aspx](http://www.nibsc.org/About_Us/Terms_and_Conditions.aspx) or upon request by the Recipient) ("Conditions") apply to the exclusion of all other terms and are hereby incorporated into this document by reference. The Recipient's attention is drawn in particular to the provisions of clause 11 of the Conditions.

#### 16. INFORMATION FOR CUSTOMS USE ONLY

<b>Country of origin for customs purposes*:</b> United Kingdom * Defined as the country where the goods have been produced and/or sufficiently processed to be classed as originating from the country of supply, for example a change of state such as freeze-drying.
<b>Net weight:</b> 0.25g
<b>Toxicity Statement:</b> Non-toxic
<b>Veterinary certificate or other statement</b> if applicable.
<b>Attached:</b> No

#### 17. CERTIFICATE OF ANALYSIS

NIBSC does not provide a Certificate of Analysis for WHO Biological Reference Materials because they are internationally recognised primary reference materials fully described in the instructions for use. The reference materials are established according to the WHO Recommendations for the preparation, characterization and establishment of international and other biological reference standards [http://www.who.int/bloodproducts/publications/TRS932Annex2\\_1nter\\_biolefstandardsrev2004.pdf](http://www.who.int/bloodproducts/publications/TRS932Annex2_1nter_biolefstandardsrev2004.pdf) (revised 2004). They are officially endorsed by the WHO Expert Committee on Biological Standardization (ECBS) based on the report of the international collaborative study which established their suitability for the intended use.