



CE Marked Material
NHS Winter multiplex Working Reagent for Nucleic Acid
Amplification Testing (NAT)
NIBSC code: 20/170-XXX
Instructions for use
(Version 5.0, Dated 12/03/2024)

This material is a self certified IVD and complies with the requirements of the "EU in vitro diagnostic medical device directive 98/79/EC"

1. INTENDED USE

This product is CE marked for use as an IVD within the UK, EU member states and EEA countries. In all other territories this product can be used for research purposes only.

The reagent is supplied to professional users, typically hospital laboratories, public health organisations and appropriate research organisations. The NIBSC NHS Winter multiplex reagent for nucleic acid amplification testing intended to be used as a run control for routine nucleic acid amplification techniques (NAT) assays for the detection of Influenza A H1N1, Influenza A H3N2, Influenza B, Respiratory Syncytial Virus A (RSV A), Respiratory Syncytial Virus B (RSV B) and the 2019 novel coronavirus SARS-CoV-2 viral contamination of clinical samples from patients who have signs and symptoms suggestive of viral respiratory disease. The control should be extracted and amplified alongside unknown samples as part of a continuing quality control programme used to monitor assay performance. The use of re-frozen or diluted product or by non-professional users may lead to inconsistent/erroneous results.

2. CAUTION

This preparation is not for administration to humans or animals in the human food chain.

The SARS-CoV-2 component of this product has been inactivated but the other components have not. For the purposes of any risk assessment this product should be categorised and handled as though it were an ACDP Hazard Group 2 pathogen.

As with all materials of biological origin, this preparation should be regarded as potentially hazardous to health. It should be used and discarded according to your own laboratory's safety procedures. Such safety procedures should include the wearing of protective gloves and avoiding the generation of aerosols. Care should be exercised in opening ampoules or vials, to avoid cuts.

3. UNITAGE

There is no unitage assigned to this control. The control should therefore be validated for use as a run control and the expected results determined by the end user for their particular NAT assay. NIBSC have determined that the Ct value of the control is approximately 30 using our in-house assay. However different extraction and amplification instruments and different assays may yield different results. Therefore it is important that each user validates this control using their own instruments and assays. Due to the slight variation between batches users are advised to re-validate their assays when using a new batch of control.

This material MUST NOT be used for any calibration purposes at all. NIBSC provide a series of dilution controls that can be used by laboratories to discover the end point sensitivity of their assay. The materials used in these dilution controls are exactly the same as used in this product 20/170.

The dilution controls are

20/180 Influenza A (H1N1) Calibrants – panel of six dilutions
20/182 Influenza A (H3N2) Calibrants – panel of six dilutions
20/184 Influenza B Calibrants – panel of six dilutions
20/186 RSV A Calibrants – panel of six dilutions
20/188 RSV B Calibrants – panel of six dilutions
20/190 SARS-CoV-2 Calibrants – panel of six dilutions

4. CONTENTS

Country of origin of biological material: United Kingdom.
Each vial is encoded 20/170-XXX and contains 1.0ml of control.
The controls consists of a whole virus preparation of
Influenza virus A H1 N1 pdm 09 (A/California/7/09)
Influenza virus A H3 N2 (A/Wyoming/3/2003)
Influenza virus B (B/Jiangsu/10/2003)
RSV A A2
RSV B 9320
SARS-CoV-2 Melbourne
diluted in a buffer comprising of PBS A and 2% foetal calf serum.

5. STORAGE

The control should be delivered in a frozen state and then be stored at or below -20°C until use. Should the material arrived in a thawed state it should be discarded and NIBSC contacted for a replacement. Material should be thawed and not refrozen. Once thawed each vial should be stored between +2°C and +8°C and then be used within five days. After this point the control should be discarded. Users are encouraged to inform NIBSC of the performance of the preparation from reviews of their data monitoring. Any user who has data supporting any deterioration in the characteristics of any reference preparation is encouraged to contact NIBSC.

6. DIRECTIONS FOR OPENING

Vials have a screw cap; an internal stopper may also be present. The cap should be removed by turning anti-clockwise. Care should be taken to prevent loss of the contents. Please note: If a stopper is present on removal of the cap, the stopper should remain in the vial or be removed with the cap.

7. USE OF MATERIAL

The control should be used directly without further dilution and extracted and amplified alongside samples under test. Best results are achieved when the entire volume of the control is extracted. It is recommended that the control be included in each assay run to monitor assay performance.

The Result Reporting System (RRS) has been developed by the National Institute for Biological Standards and Control (NIBSC) for the data monitoring of its serology and NAT quality control (QC) reagents. These include the Quality Control Reagent Unit (QCRU) and Clinical Virology Network (CVN) reagents. The system has been successfully running for serology assays for several years, collecting thousands of data points a year; and more recently for data derived from reagents used in Nucleic Acid-based Technologies (NAT) assays. Users are encouraged to sign up here: <https://nibsc.org/products/rrs>

8. STABILITY

Reference materials are held at NIBSC within assured, temperature-controlled storage facilities. Reference Materials should be stored on receipt as indicated on the label.
NIBSC follows the policy of WHO with respect to its reference materials.



9. REFERENCES

10. ACKNOWLEDGEMENTS

1. This product has been made at NIBSC in collaboration with Dr Mel Smith and Dr Angela Douglas and the Clinical Virology network. NIBSC acknowledge the assistance of those laboratories that provided initial data for this product especially those scientists working in PHE and NHS labs in Bristol, Birmingham, Cambridge, Ipswich, Newcastle, Norwich and Nottingham.

2. EC REP: Advena Ltd. Tower Business Centre, 2nd Flr., Tower Street, Swatar, BKR 4013 Malta

11. FURTHER INFORMATION

Further information can be obtained as follows;

This material: enquiries@nibsc.org

WHO Biological Standards:

<http://www.who.int/biologicals/en/>

JCTLM Higher order reference materials:

<http://www.bipm.org/en/committees/jc/jctlm/>

Derivation of International Units:

http://www.nibsc.org/standardisation/international_standards.aspx

Ordering standards from NIBSC:

<http://www.nibsc.org/products/ordering.aspx>

NIBSC Terms & Conditions:

http://www.nibsc.org/terms_and_conditions.aspx

12. CUSTOMER FEEDBACK

Customers are encouraged to provide feedback on the suitability or use of the material provided or other aspects of our service. Please send any comments to enquiries@nibsc.org

13. CITATION

In all publications, including data sheets, in which this material is referenced, it is important that the preparation's title, its status, the NIBSC code number, and the name and address of NIBSC are cited and cited correctly.

14. MATERIAL SAFETY SHEET

Classification in accordance with Directive 2000/54/EC, Regulation (EC) No 1272/2008: Not applicable or not classified

Physical and Chemical properties	
Physical appearance: Frozen Liquid	Corrosive: No
Stable: Yes	Oxidising: No
Hygroscopic: No	Irritant: No
Flammable: No	Handling: See caution, Section 2
Other (specify): N/A	
Toxicological properties	
Effects of inhalation:	Not established, avoid inhalation
Effects of ingestion:	Not established, avoid ingestion
Effects of skin absorption:	Not established, avoid contact with skin
Suggested First Aid	
Inhalation:	Seek medical advice
Ingestion:	Seek medical advice

Contact with eyes:	Wash with copious amounts of water. Seek medical advice
Contact with skin:	Wash thoroughly with water.
Action on Spillage and Method of Disposal	
Spillage of ampoule contents should be taken up with absorbent material wetted with an appropriate disinfectant. Rinse area with an appropriate disinfectant followed by water. Absorbent materials used to treat spillage should be treated as biological waste.	

15. LIABILITY AND LOSS

In the event that this document is translated into another language, the English language version shall prevail in the event of any inconsistencies between the documents.

Unless expressly stated otherwise by NIBSC, NIBSC's Standard Terms and Conditions for the Supply of Materials (available at http://www.nibsc.org/About_Us/Terms_and_Conditions.aspx or upon request by the Recipient) ("Conditions") apply to the exclusion of all other terms and are hereby incorporated into this document by reference. The Recipient's attention is drawn in particular to the provisions of clause 11 of the Conditions.

16. INFORMATION FOR CUSTOMS USE ONLY

Country of origin for customs purposes*: United Kingdom
* Defined as the country where the goods have been produced and/or sufficiently processed to be classed as originating from the country of supply, for example a change of state such as freeze-drying.
Net weight: 1g
Toxicity Statement: Non-toxic
Veterinary certificate or other statement if applicable.
Attached: No