



**WHO International Standard**  
**First WHO International Standard for Plasmodium falciparum**  
**antigens**  
**NIBSC code: 16/376**  
**Instructions for use**  
**(Version 1.0, Dated 22/11/2017)**

**1. INTENDED USE**

The 1st WHO International Standard for Plasmodium falciparum antigens consists of a freeze-dried preparation of culture-derived P. falciparum parasites of the W2 strain. The standard has been lyophilised in 0.5 mL aliquots and stored at -70°C. This preparation has been assessed in a collaborative study for its suitability for use in a range of P. falciparum antigen detection tests. The collaborative study report contains full details of the reactivity of 16/376 [1]. The intended uses of this International Standard are to allow the standardisation, and evaluation of performance and sensitivity, of P. falciparum HRP2 and pLDH antigen detection tests.

**2. CAUTION**

**This preparation is not for administration to humans or animals in the human food chain.**

The preparation contains material of human origin, and either the final product or the source materials, from which it is derived, have been tested and found negative for HBsAg, anti-HIV and HCV RNA. As with all materials of biological origin, this preparation should be regarded as potentially hazardous to health. It should be used and discarded according to your own laboratory's safety procedures. Such safety procedures should include the wearing of protective gloves and avoiding the generation of aerosols. Care should be exercised in opening ampoules or vials, to avoid cuts.

**3. UNITAGE**

This material has been assigned a unitage of 1000 International Units of HRP2 per ampoule and 1000 International Units of pLDH per ampoule.

**4. CONTENTS**

Country of origin of biological material: United States of America.  
Parasite extract prepared from a P. falciparum in vitro culture of the W2 strain was diluted (1:75) with 10 mM Tris pH 7.4, 5% trehalose, 1 mM EDTA solution and filled at 0.5 mL/ampoule.

**5. STORAGE**

This preparation should be stored at -20°C or below on receipt. If the material is to be stored for more than 3 months prior to use storage at -70°C is recommended.

**Please note: because of the inherent stability of lyophilized material, NIBSC may ship these materials at ambient temperature.**

**6. DIRECTIONS FOR OPENING**

DIN ampoules have an 'easy-open' coloured stress point, where the narrow ampoule stem joins the wider ampoule body. Various types of ampoule breaker are available commercially. To open the ampoule, tap the ampoule gently to collect material at the bottom (labelled) end and follow manufactures instructions provided with the ampoule breaker.

**7. USE OF MATERIAL**

**No attempt should be made to weigh out any portion of the freeze-dried material prior to reconstitution**

This material is supplied lyophilised and before use should be reconstituted in 0.5 mL of whole blood. Reconstituted material should be used on the day of reconstitution.

**8. STABILITY**

Reference materials are held at NIBSC within assured, temperature-controlled storage facilities. Reference Materials should be stored on receipt as indicated on the label.

NIBSC follows the policy of WHO with respect to its reference materials.

**9. REFERENCES**

1. Lynne M. Harris, Ana Campillo, Peter Rigsby, Eleanor Atkinson, Ryaka Poonawala, Michael Aidoo, John Saldanha, Iveth J. Gonzalez, Paul W. Bowyer, and the Collaborative Study Group (2017). Collaborative study to evaluate the proposed First World Health Organization International Standard for Plasmodium falciparum antigens. WHO/BS/2017.2329.

**10. ACKNOWLEDGEMENTS**

This material was developed in collaboration with the Foundation for Innovative New Diagnostics (FIND). We gratefully acknowledge the significant contributions made by the participating laboratories and scientists involved in the collaborative study. We would like to thank the US CDC for providing the parasite source material for this International Standard.

**11. FURTHER INFORMATION**

Further information can be obtained as follows;

This material: [enquiries@nibsc.org](mailto:enquiries@nibsc.org)

WHO Biological Standards:

<http://www.who.int/biologicals/en/>

JCTLM Higher order reference materials:

<http://www.bipm.org/en/committees/jc/jctlm/>

Derivation of International Units:

[http://www.nibsc.org/standardisation/international\\_standards.aspx](http://www.nibsc.org/standardisation/international_standards.aspx)

Ordering standards from NIBSC:

<http://www.nibsc.org/products/ordering.aspx>

NIBSC Terms & Conditions:

[http://www.nibsc.org/terms\\_and\\_conditions.aspx](http://www.nibsc.org/terms_and_conditions.aspx)

**12. CUSTOMER FEEDBACK**

Customers are encouraged to provide feedback on the suitability or use of the material provided or other aspects of our service. Please send any comments to [enquiries@nibsc.org](mailto:enquiries@nibsc.org)

**13. CITATION**

In all publications, including data sheets, in which this material is referenced, it is important that the preparation's title, its status, the NIBSC code number, and the name and address of NIBSC are cited and cited correctly.

**14. MATERIAL SAFETY SHEET**

Classification in accordance with Directive 2000/54/EC, Regulation (EC) No 1272/2008: Not applicable or not classified

| Physical and Chemical properties  |  |
|-----------------------------------|--|
| Physical appearance: Freeze-dried | Corrosive: No                            |
| Stable: Yes                       | Oxidising: No                            |
| Hygroscopic: No                   | Irritant: No                             |
| Flammable: No                     | Handling: See caution, Section 2         |
| Other (specify):                  | Contains material of human origin        |
| Toxicological properties          |  |
| Effects of inhalation:            | Not established, avoid inhalation        |
| Effects of ingestion:             | Not established, avoid ingestion         |
| Effects of skin absorption:       | Not established, avoid contact with skin |



| <b>Suggested First Aid</b>   |   |
|--|---|
| Inhalation:  | Seek medical advice                                     |
| Ingestion:   | Seek medical advice                                     |
| Contact with eyes:   | Wash with copious amounts of water. Seek medical advice |
| Contact with skin:   | Wash thoroughly with water.                             |
| <b>Action on Spillage and Method of Disposal</b>   |   |
| Spillage of ampoule contents should be taken up with absorbent material wetted with an appropriate disinfectant. Rinse area with an appropriate disinfectant followed by water.<br>Absorbent materials used to treat spillage should be treated as biological waste. |   |

#### 15. LIABILITY AND LOSS

In the event that this document is translated into another language, the English language version shall prevail in the event of any inconsistencies between the documents.

Unless expressly stated otherwise by NIBSC, NIBSC's Standard Terms and Conditions for the Supply of Materials (available at [http://www.nibsc.org/About\\_Us/Terms\\_and\\_Conditions.aspx](http://www.nibsc.org/About_Us/Terms_and_Conditions.aspx) or upon request by the Recipient) ("Conditions") apply to the exclusion of all other terms and are hereby incorporated into this document by reference. The Recipient's attention is drawn in particular to the provisions of clause 11 of the Conditions.

#### 16. INFORMATION FOR CUSTOMS USE ONLY

|   |
|---|
| <b>Country of origin for customs purposes*:</b> United Kingdom<br>* Defined as the country where the goods have been produced and/or sufficiently processed to be classed as originating from the country of supply, for example a change of state such as freeze-drying. |
| <b>Net weight:</b> 0.5211g  |
| <b>Toxicity Statement:</b> Non-toxic  |
| <b>Veterinary certificate or other statement</b> if applicable.   |
| <b>Attached:</b> No   |

#### 17. CERTIFICATE OF ANALYSIS

NIBSC does not provide a Certificate of Analysis for WHO Biological Reference Materials because they are internationally recognised primary reference materials fully described in the instructions for use. The reference materials are established according to the WHO Recommendations for the preparation, characterization and establishment of international and other biological reference standards [http://www.who.int/bloodproducts/publications/TRS932Annex2\\_Inter\\_bi\\_olefstandardsrev2004.pdf](http://www.who.int/bloodproducts/publications/TRS932Annex2_Inter_bi_olefstandardsrev2004.pdf) (revised 2004). They are officially endorsed by the WHO Expert Committee on Biological Standardization (ECBS) based on the report of the international collaborative study which established their suitability for the intended use.