



**WHO International Standard  
1st IS for Mycobacterium tuberculosis (H37Rv) DNA for NAT-  
based assays  
NIBSC code: 20/152  
Instructions for use  
(Version 3.0, Dated 05/11/2021)**

**1. INTENDED USE**

This inactivated Mycobacterium tuberculosis strain H37Rv preparation was generated by NIBSC, U.K. in 2020. The strain H37Rv was sourced from National Collection of Type Cultures in Public Health England, U.K. in 2018. It was established as the 1<sup>st</sup> WHO International Standard for Mycobacterium tuberculosis (H37Rv) DNA for NAT-based assays in 2021. The intended uses of this material are for calibration of secondary or in-house reference materials used in the assays for the molecular detection of M. tuberculosis DNA. It may also be used for assay validation and monitoring the limit of detection of rapid diagnostic tests.

**2. CAUTION**

**This preparation is not for administration to humans or animals in the human food chain.**

The material is not of human or bovine origin. As with all materials of biological origin, this preparation should be regarded as potentially hazardous to health. It should be used and discarded according to your own laboratory's safety procedures. Such safety procedures should include the wearing of protective gloves and avoiding the generation of aerosols. Care should be exercised in opening ampoules or vials, to avoid cuts.

**3. UNITAGE**

This preparation contains an arbitrary unitage of 6.3 log<sub>10</sub> (or 2 million) IU per vial.

**4. CONTENTS**

Country of origin of biological material: United Kingdom.  
Each vial contains the residue after freeze-drying of 0.5 ml of inactivated M. tuberculosis (H37Rv) suspension in TE buffer (10 mM Tris-HCl pH 7.4, 1 mM EDTA) containing 2% D-(+)-Trehalose dehydrate and 0.2% Tween-80. There are approximately 6.3 log<sub>10</sub> (or 2 million) genome copies with a standard deviation of 0.7 per vial as estimated from a collaborative study (see reference in section 9).

**5. STORAGE**

This preparation should be stored at -20°C for long-term storage.  
**Please note: because of the inherent stability of lyophilized material, NIBSC may ship these materials at ambient temperature.**

**6. DIRECTIONS FOR OPENING**

Vials have a screw cap; an internal stopper may also be present. The cap should be removed by turning anti-clockwise. Care should be taken to prevent loss of the contents. Please note: If a stopper is present on removal of the cap, the stopper should remain in the vial or be removed with the cap.

**7. USE OF MATERIAL**

**No attempt should be made to weigh out any portion of the freeze-dried material prior to reconstitution**

The content of the vial can be reconstituted in sterile deionised water or appropriate buffer. The reconstituted preparation can be stored at 4°C for up to 1 week or at -20°C for up to 4 weeks.

**8. STABILITY**

Reference materials are held at NIBSC within assured, temperature-controlled storage facilities. Reference Materials should be stored on receipt as indicated on the label.

NIBSC follows the policy of WHO with respect to its reference materials.

**9. REFERENCES**

WHO/BS/2021.2403

**10. ACKNOWLEDGEMENTS**

We gratefully acknowledge the important contributions of the collaborative study participants.

**11. FURTHER INFORMATION**

Further information can be obtained as follows;  
This material: [enquiries@nibsc.org](mailto:enquiries@nibsc.org)  
WHO Biological Standards:  
<http://www.who.int/biologicals/en/>  
JCTLM Higher order reference materials:  
<http://www.bipm.org/en/committees/jc/jctlm/>  
Derivation of International Units:  
[http://www.nibsc.org/standardisation/international\\_standards.aspx](http://www.nibsc.org/standardisation/international_standards.aspx)  
Ordering standards from NIBSC:  
<http://www.nibsc.org/products/ordering.aspx>  
NIBSC Terms & Conditions:  
[http://www.nibsc.org/terms\\_and\\_conditions.aspx](http://www.nibsc.org/terms_and_conditions.aspx)

**12. CUSTOMER FEEDBACK**

Customers are encouraged to provide feedback on the suitability or use of the material provided or other aspects of our service. Please send any comments to [enquiries@nibsc.org](mailto:enquiries@nibsc.org)

**13. CITATION**

In all publications, including data sheets, in which this material is referenced, it is important that the preparation's title, its status, the NIBSC code number, and the name and address of NIBSC are cited and cited correctly.

**14. MATERIAL SAFETY SHEET**

Classification in accordance with Directive 2000/54/EC, Regulation (EC) No 1272/2008: Not applicable or not classified

Physical and Chemical properties	
Physical appearance: Freeze dried white cake	Corrosive: No
Stable: Yes	Oxidising: No
Hygroscopic: Yes	Irritant: No
Flammable: No	Handling: See caution, Section 2
Other (specify):	Contains heat inactivated Mycobacterium tuberculosis strain H37Rv as non-infectious material
Toxicological properties	
Effects of inhalation:	Not established, avoid inhalation
Effects of ingestion:	Not established, avoid ingestion
Effects of skin absorption:	Not established, avoid contact with skin
Suggested First Aid	
Inhalation:	Seek medical advice
Ingestion:	Seek medical advice
Contact with eyes:	Wash with copious amounts of water. Seek medical advice
Contact with skin:	Wash thoroughly with water.



**Action on Spillage and Method of Disposal**

Spillage of vial contents should be taken up with absorbent material wetted with an appropriate disinfectant. Rinse area with an appropriate disinfectant followed by water.  
Absorbent materials used to treat spillage should be treated as biological waste.

**15. LIABILITY AND LOSS**

In the event that this document is translated into another language, the English language version shall prevail in the event of any inconsistencies between the documents.

Unless expressly stated otherwise by NIBSC, NIBSC's Standard Terms and Conditions for the Supply of Materials (available at [http://www.nibsc.org/About\\_Us/Terms\\_and\\_Conditions.aspx](http://www.nibsc.org/About_Us/Terms_and_Conditions.aspx) or upon request by the Recipient) ("Conditions") apply to the exclusion of all other terms and are hereby incorporated into this document by reference. The Recipient's attention is drawn in particular to the provisions of clause 11 of the Conditions.

**16. INFORMATION FOR CUSTOMS USE ONLY**

**Country of origin for customs purposes\*:** United Kingdom

\* Defined as the country where the goods have been produced and/or sufficiently processed to be classed as originating from the country of supply, for example a change of state such as freeze-drying.

**Net weight:** < 10 mg

**Toxicity Statement:** Non-toxic

**Veterinary certificate or other statement** if applicable.

**Attached:** No

**17. CERTIFICATE OF ANALYSIS**

NIBSC does not provide a Certificate of Analysis for WHO Biological Reference Materials because they are internationally recognised primary reference materials fully described in the instructions for use. The reference materials are established according to the WHO Recommendations for the preparation, characterization and establishment of international and other biological reference standards [http://www.who.int/bloodproducts/publications/TRS932Annex2\\_Inter\\_biologicalstandardsrev2004.pdf](http://www.who.int/bloodproducts/publications/TRS932Annex2_Inter_biologicalstandardsrev2004.pdf) (revised 2004). They are officially endorsed by the WHO Expert Committee on Biological Standardization (ECBS) based on the report of the international collaborative study which established their suitability for the intended use.