

Non WHO Reference Material
Anti-Meningococcal Serotype P3.15 Monoclonal Antibody
NIBSC code: 02/310
Instructions for use
(Version 5.0, Dated 06/11/2013)

This material is not for in vitro diagnostic use.

1. INTENDED USE

For use as a typing reagent.

2. CAUTION

This preparation is not for administration to humans or animals in the human food chain.

The preparation contains material of bovine origin that is certified to be obtained from animals taken from a closed herd in the female line since 1980, in which no animal has been clinically suspected of having BSE and which has not been fed rations containing ruminant derived protein during that period. As with all materials of biological origin, this preparation should be regarded as potentially hazardous to health. It should be used and discarded according to your own laboratory's safety procedures. Such safety procedures should include the wearing of protective gloves and avoiding the generation of aerosols. Care should be exercised in opening ampoules or vials, to avoid cuts.

UNITAGE

Units of activity have not been assigned to this material. Please refer to table on page 2 for recommended working concentrations.

4. CONTENTS

Country of origin of biological material: United Kingdom.

Each ampoule contains the freeze-dried powder from 1ml of cell culture supernatant concentrated approximately 15 fold. Antibody is of murine origin.

5. STORAGE

Store freeze dried ampoules and reconstituted aliquots at -20°C.

Please note: because of the inherent stability of lyophilized material, NIBSC may ship these materials at ambient temperature.

6. DIRECTIONS FOR OPENING

DIN ampoules have an 'easy-open' coloured stress point, where the narrow ampoule stem joins the wider ampoule body. Various types of ampoule breaker are available commercially. To open the ampoule, tap the ampoule gently to collect material at the bottom (labelled) end and follow manufactures instructions provided with the ampoule breaker.

7. USE OF MATERIAL

No attempt should be made to weigh out any portion of the freeze-dried material prior to reconstitution

Resuspend each din ampoule with 1 ml distilled water. Ensure the entire content of each ampoule is fully resuspended.

8. STABILITY

Reference materials are held at NIBSC within assured, temperature-controlled storage facilities. Reference Materials should be stored on receipt as indicated on the label.

The recommended working concentrations were correct at the time of manufacture – no information is available on long term stability. Stability of the reconstituted material should be determined by the user. Users who have data supporting any deterioration in the characteristics of any reference preparation are encouraged to contact NIBSC.

9. REFERENCES

Zollinger, W.D., E.E. Moran, H. Connelly, R.E. Mandrell and B. Brandt, 1984. Monoclonal antibodies to serotype 2 and serotype 15 outer membrane proteins of *Neisseria meningitidis* and their uses in serotyping. Infection and Immunity 46:260-266.

10. ACKNOWLEDGEMENTS

This material was produced from the hybridoma cell line 8B5-5-G9. Provided by Dr W.D. Zollinger of the Walter Reed Army Institute of Research, Washington D.C., U.S.A..

11. FURTHER INFORMATION

Further information can be obtained as follows;

This material: enquiries@nibsc.org

WHO Biological Standards:

http://www.who.int/biologicals/en/

JCTLM Higher order reference materials:

http://www.bipm.org/en/committees/jc/jctlm/

Derivation of International Units:

http://www.nibsc.org/standardisation/international_standards.aspx

Ordering standards from NIBSC:

http://www.nibsc.org/products/ordering.aspx

NIBSC Terms & Conditions:

http://www.nibsc.org/terms_and_conditions.aspx

12. CUSTOMER FEEDBACK

Customers are encouraged to provide feedback on the suitability or use of the material provided or other aspects of our service. Please send any comments to enquiries@nibsc.org

13. CITATION

In all publications, including data sheets, in which this material is referenced, it is important that the preparation's title, its status, the NIBSC code number, and the name and address of NIBSC are cited and cited correctly.

14. MATERIAL SAFETY SHEET

Classification in accordance with Directive 2000/54/EC, Regulation (EC) No 1272/2008: Not applicable or not classified

110 1212/2000. Not applicable of flot diagomed								
Physical and Chemical properties								
Physical	Corrosive:	No						
appearance:								
Freeze								
dried powder								
Stable:	Oxidising:	No						
Yes								
Hygroscopic:	Irritant:	No						
No								
Flammable:	Handling:	See caution, Section 2						
No								
Other (specify): N/A								
Toxicological properties								
Toxicological properties								
Effects of inhalation: Not established, avoid inhalation								
Effects of ingestion: Not established, avoid ingestion								
Effects of skin absorption: Not established,								
avoid contact with skin								
Suggested First Aid								
ouggested i list Aid								
Inhalation: Seek medical advice								
Ingestion: Seek medical advice								
Contact with eyes: Wash with copious amounts of water								
Seek medical advice								

National Institute for Biological Standards and Control,

Potters Bar, Hertfordshire, EN6 3QG. T +44 (0)1707 641000, nibsc.org WHO International Laboratory for Biological Standards, UK Official Medicines Control Laboratory



Contact with skin: Wash thoroughly with water.

Action on Spillage and Method of Disposal

Spillage of ampoule contents should be taken up with absorbent material wetted with an appropriate disinfectant. Rinse area with an appropriate disinfectant followed by water.

Absorbent materials used to treat spillage should be treated as biological waste.

15. LIABILITY AND LOSS

In the event that this document is translated into another language, the English language version shall prevail in the event of any inconsistencies between the documents.

Unless expressly stated otherwise by NIBSC, NIBSC's Standard Terms and Conditions for the Supply of Materials (available at http://www.nibsc.org/About_Us/Terms_and_Conditions.aspx or upon

request by the Recipient) ("Conditions") apply to the exclusion of all other terms and are hereby incorporated into this document by reference. The Recipient's attention is drawn in particular to the provisions of clause 11 of the Conditions.

16. INFORMATION FOR CUSTOMS USE ONLY

Country of origin for customs purposes*: United Kingdom

* Defined as the country where the goods have been produced and/or sufficiently processed to be classed as originating from the country of supply, for example a change of state such as freeze-drying.

Net weight: 4.5g.

Toxicity Statement: Toxicity not assessed

Veterinary certificate or other statement if applicable.

Attached: No

3. UNITAGE continued

The material has been tested for use in whole cell dot-blot and whole cell ELISA as follows:

Specificity	Source of mAb ¹	NIBSC hybridoma stock number ²	Isotype	Resuspension	Concentration of reconstituted stock to use in dot-blot	mAb 'name'
Serosubtype P3.15	Zollinger 8B5-5-G9	4061	IgG2a	Each vial should be resuspended with 1ml distilled water	1 in 250	2-1-P15

¹ Source of mAb: This indicates the laboratory where the hybridoma cells were isolated and the hybridoma clone designation (this is also thename used in the literature for the NVI hybridomas).

National Institute for Biological Standards and Control,

Potters Bar, Hertfordshire, EN6 3QG. T +44 (0)1707 641000, nibsc.org WHO International Laboratory for Biological Standards, UK Official Medicines Control Laboratory

NIBSC hybridoma stock number: this number was assigned at NIBSC when the hybridoma cells were received and refers to cells rather than antibody