

#### WHO International Standard Anti-D Immunoglobulin NIBSC code: 16/332 Instructions for use (Version 3.0, Dated 25/07/2019)

### 1. INTENDED USE

Preparation 16/332, the 3<sup>rd</sup> International Standard for anti-D immunoglobulin, is intended to be used as a standard in potency assays of anti-D immunoglobulin products.

# 2. CAUTION

# This preparation is not for administration to humans or animals in the human food chain.

The preparation contains material of human origin, and either the final product or the source materials, from which it is derived, have been tested and found negative for HBsAg, anti-HIV and HCV RNA. As with all materials of biological origin, this preparation should be regarded as potentially hazardous to health. It should be used and discarded according to your own laboratory's safety procedures. Such safety procedures should include the wearing of protective gloves and avoiding the generation of aerosols. Care should be exercised in opening ampoules or vials, to avoid cuts.

# 3. UNITAGE

EACH AMPOULE CONTAINS 297 INTERNATIONAL UNITS OF ANTI-D IMMUNOGLOBULIN.

# 4. CONTENTS

Country of origin of biological material: UK.

Preparation 16/332 consists of anti-D immunoglobulin product manufactured from USA plasma. Each ampoule contains the lyophilized residue of 1 ml of 2% (w/v) IgG in 0.25M glycine, 0.4% NaCl, pH 6.5.

#### 5. STORAGE

Store unopened ampoules at -20°C or below.

Please note: because of the inherent stability of lyophilized material, NIBSC may ship these materials at ambient temperature.

### 6. DIRECTIONS FOR OPENING

DIN ampoules have an 'easy-open' coloured stress point, where the narrow ampoule stem joins the wider ampoule body. Various types of ampoule breaker are available commercially. To open the ampoule, tap the ampoule gently to collect material at the bottom (labelled) end and follow manufactures instructions provided with the ampoule breaker.

#### 7. USE OF MATERIAL

No attempt should be made to weigh out any portion of the freeze-dried material prior to reconstitution

Reconstitute the ampoule contents with 1.0 ml distilled or deionised water containing 0.02% sodium azide.

#### 8. STABILITY

It is the policy of WHO not to assign an expiry date to its international reference materials. They remain valid with the assigned potency and status until withdrawn or amended.

Reference materials are held at NIBSC within assured, temperaturecontrolled storage facilities. Reference Materials should be stored on receipt as indicated on the label. For information specific to a particular biological standard, contact NIBSC.

Accelerated degradation studies on 16/332 indicate that the lyophilised material will be adequately stable at -20°C with respect to anti-D

National Institute for Biological Standards and Control,

Potters Bar, Hertfordshire, EN6 3QG. T +44 (0)1707 641000, nibsc.org WHO International Laboratory for Biological Standards, UK Official Medicines Control Laboratory

potency. Once reconstituted, diluted or aliquoted, users should determine the stability of the material according to their own method of preparation, storage and use. However, tests at NIBSC have shown that reconstituted 16/332 is stable for at least 1 month at 4°C in a tightly-lidded tube in the presence of 0.02% (w/v) sodium azide.

NIBSC follows the policy of WHO with respect to its reference materials.

Users who have data supporting any deterioration in the characteristics of any reference preparation are encouraged to contact NIBSC.

### 9. REFERENCES

Fox, B. , Sharp, G. , Atkinson, E. , Roberts, G. , Rigsby, P. and Studholme, L. (2019), The third international standard for anti-D immunoglobulin: international collaborative study to evaluate candidate preparations. Vox Sang. https://doi.org/10.1111/vox.12822

# 10. ACKNOWLEDGEMENTS

We thank the participants of the collaborative study to assign unitage.

# 11. FURTHER INFORMATION

Further information can be obtained as follows; This material: enquiries@nibsc.org WHO Biological Standards: http://www.who.int/biologicals/en/ JCTLM Higher order reference materials: http://www.bipm.org/en/committees/jc/jctlm/ Derivation of International Units: http://www.nibsc.org/standardisation/international\_standards.aspx Ordering standards from NIBSC: http://www.nibsc.org/products/ordering.aspx NIBSC Terms & Conditions: http://www.nibsc.org/terms\_and\_conditions.aspx

# 12. CUSTOMER FEEDBACK

Customers are encouraged to provide feedback on the suitability or use of the material provided or other aspects of our service. Please send any comments to enquiries@nibsc.org

#### 13. CITATION

In all publications, including data sheets, in which this material is referenced, it is important that the preparation's title, its status, the NIBSC code number, and the name and address of NIBSC are cited and cited correctly.

#### 14. MATERIAL SAFETY SHEET

Classification in accordance with Directive 2000/54/EC, Regulation (EC) No 1272/2008: Not applicable or not classified

Physical and Chemical properties					
Physical appearance:		Corrosive:	No		
Lyophilisate					
Stable:	Yes		Oxidising:	No	
Hygroscopic:	No		Irritant:	No	
Flammable:	No		Handling: Se	e caution, Section 2	
Other (specify): Contains material of human origin					
Toxicological properties					
Effects of inhalation:		Not established, avoid inhalation			
Effects of ingestion:		Not established, avoid ingestion			
Effects of skin absorption:		Not established, avoid contact with skin			
Suggested First Aid					
Inhalation: Seek medical advice					
Ingestion: Seek medical advice					



# Medicines & Healthcare products Regulatory Agency



Action on Spillage and Method of Disposal					
Contact with skin:	Wash thoroughly with water.				
Contact with eyes:	Wash with copious amounts of water. Seek medical advice				

#### Action on Spillage and Method of Disposal

Spillage of contents should be taken up with absorbent material wetted with an appropriate disinfectant. Rinse area with an appropriate disinfectant followed by water.

Absorbent materials used to treat spillage should be treated as biological waste.

# 15. LIABILITY AND LOSS

In the event that this document is translated into another language, the English language version shall prevail in the event of any inconsistencies between the documents.

Unless expressly stated otherwise by NIBSC, NIBSC's Standard Terms and Conditions for the Supply of Materials (available at http://www.nibsc.org/About\_Us/Terms\_and\_Conditions.aspx or upon request by the Recipient) ("Conditions") apply to the exclusion of all other terms and are hereby incorporated into this document by reference. The Recipient's attention is drawn in particular to the provisions of clause 11 of the Conditions.

# 16. INFORMATION FOR CUSTOMS USE ONLY

Country of origin for customs purposes\*: United Kingdom \* Defined as the country where the goods have been produced and/or sufficiently processed to be classed as originating from the country of supply, for example a change of state such as freeze-drying. Net weight: 0.04g Toxicity Statement: Toxicity not assessed Veterinary certificate or other statement if applicable.

Attached: No

# 17. CERTIFICATE OF ANALYSIS

NIBSC does not provide a Certificate of Analysis for WHO Biological Reference Materials because they are internationally recognised primary reference materials fully described in the instructions for use. The reference materials are established according to the WHO Recommendations for the preparation, characterization and establishment of international and other biological reference standards http://www.who.int/bloodproducts/publications/TRS932Annex2\_Inter\_bi olefstandardsrev2004.pdf (revised 2004). They are officially endorsed by the WHO Expert Committee on Biological Standardization (ECBS) based on the report of the international collaborative study which established their suitability for the intended use.

National Institute for Biological Standards and Control, Potters Bar, Hertfordshire, EN6 3QG. T +44 (0)1707 641000, nibsc.org WHO International Laboratory for Biological Standards, UK Official Medicines Control Laboratory

