

1. INTENDED USE

This consists of a batch of ampoules (coded 75/569) containing purified alpha subunit of human chorionic gonadotrophin (hCG-alpha). As authorized by the WHO Expert Committee on Biological Standardization (WHO ECBS, 1975), this preparation was established by the National Institute for Biological Standards and Control as the International Reference Preparation of Chorionic Gonadotrophin Alpha Subunit, Human, for Immunoassay and assigned a potency of 70 international units per ampoule (on the basis that one International Unit is the activity contained in one microgram of the pure subunit preparation).

For further details of this Standard and of its collaborative study see Storring et al (1980). For details of the purification and characterization of the hCG-alpha used to prepare the Standard see Canfield and Ross (1976).

2. CAUTION

This preparation is not for administration to humans or animals in the human food chain.

The preparation contains material of human origin, and either the final product or the source materials, from which it is derived, have been tested and found negative for HBsAg, anti-HIV and HCV RNA.

However, as with all preparations of human origin, this material cannot be assumed to be free from infectious agents. Suitable precautions should be taken in the use and disposal of the ampoule and its contents. Such safety procedures probably will include the wearing of protective gloves and avoiding the generation of aerosols. Care should be exercised in opening ampoules or vials, to avoid cuts.

3. UNITAGE

Each ampoule contains 70 INTERNATIONAL UNITS (by definition).

4. CONTENTS

Country of origin of biological material: United Kingdom.

Each ampoule contains the freeze-dried residue of 0.5 ml of a solution which contained:

HCG-Alpha subunit approx 70µg
Human plasma albumin " 5mg
sodium chloride " 445µg
acetic acid " 300µg

Nitrogen gas at slightly less than atmospheric pressure.

5. STORAGE


Unopened ampoules should be stored at -20°C.

Please note: because of the inherent stability of lyophilized material, NIBSC may ship these materials at ambient temperature.

6. DIRECTIONS FOR OPENING

Tap the ampoule gently to collect the material at the bottom (labelled) end. Ensure ampoule is scored all round at the narrow part of the neck, with a diamond or tungsten carbide tipped glass knife file or other suitable implement before attempting to open. Place the ampoule in the ampoule opener, positioning the score at position 'A'; shown in the diagram below. Surround the ampoule with cloth or layers of tissue paper. Grip the ampoule and holder in the hand and squeeze at point 'B'. The ampoule will snap open. Take care to avoid cuts and projectile glass fragments that enter eyes. Take care that no material is lost from the ampoule and that no glass falls into the ampoule.





Side view of ampoule opening device containing an ampoule positioned ready to open. 'A' is the score mark and 'B' the point of applied pressure.

7. USE OF MATERIAL

For practical purposes each ampoule contains the same amount of the same materials. Dissolve all the contents in a known amount of buffer solution. No attempt should be made to weigh portions of the freeze-dried powder.

For economy of use the solution can be kept for several months if an anti-bacterial preservative is added and the solution is subdivided into several small containers, which are frozen rapidly to below -70°C and then stored below -30°C in the dark; repeated freezing and thawing should be avoided. If extensive dilutions are prepared, a carrier protein (0.1% w/v) should be added, which is free of peptidase.

The material has not been sterilized and contains no bacteriostat.

8. STABILITY

NIBSC follows the policy of WHO with respect to its reference materials. It is the policy of WHO not to assign an expiry date to their international reference materials. They remain valid with the assigned potency and status until withdrawn or amended. Reference materials are held at NIBSC within assured, temperature-controlled storage facilities. Unopened ampoules should be stored on receipt as indicated on the label. In addition, once reconstituted, diluted or aliquoted, users should determine the stability of the material according to their own method of preparation, storage and use.

9. REFERENCES

Canfield, R.E. and G.T. Ross (1976). Bull Wld Hlth Org., **54**, 463-470. Storring P L, Gaines Das R E & Bangham D R (1980). J. Endocrinol. **84**, 295-310.

WHO Expert Committee on Biological Standardization (1975). Wld Hlth Org. Tech. Rep. Ser. No. 565.

10. ACKNOWLEDGEMENTS

The Center for Population Research, U.S.A. and the Reproductive Research Branch of the National Institute of Child Health and Development, U.S.A. for providing the hCG- α ; Drs R E Canfield and G T Ross and their colleagues for the purification and characterization of the hCG- α ; Dr R M Lequin for immunoassays; and Dr P J Campbell for ampouling.

11. FURTHER INFORMATION

Further information can be obtained as follows;

This material: enquiries@nibsc.org

WHO Biological Standards:

http://www.who.int/biologicals/en/

JCTLM Higher order reference materials:

http://www.bipm.org/en/committees/jc/jctlm/

Derivation of International Units:

http://www.nibsc.org/standardisation/international_standards.aspx

Ordering standards from NIBSC:

http://www.nibsc.org/products/ordering.aspx

NIBSC Terms & Conditions:

http://www.nibsc.org/terms_and_conditions.aspx

12. CUSTOMER FEEDBACK

Customers are encouraged to provide feedback on the suitability or use of the material provided or other aspects of our service. Please send any comments to enquiries@nibsc.org

Potters Bar, Hertfordshire, EN6 3QG. T +44 (0)1707 641000, nibsc.org WHO International Laboratory for Biological Standards, UK Official Medicines Control Laboratory

NIBSC Confidence in Biological Medicines

13. CITATION

In all publications, including data sheets, in which this material is referenced, it is important that the preparation's title, its status, the NIBSC code number, and the name and address of NIBSC are cited and cited correctly.

16. INFORMATION FOR CUSTOMS USE ONLY

Country of origin for customs purposes*: United Kingdom

* Defined as the country where the goods have been produced and/or sufficiently processed to be classed as originating from the country of supply, for example a change of state such as freeze-drying.

Net weight: 6mg

Toxicity Statement: Non-toxic

Veterinary certificate or other statement if applicable.

Attached: No

14. MATERIAL SAFETY SHEET

Classification in accordance with Directive 2000/54/EC, Regulation (EC) No 1272/2008: Not applicable or not classified

(Le) No 12/2/2000. Not applicable of not classified		
Physical and Chemical properties		
Physical	Corrosive:	No
appearance: Solid		
Stable:	Oxidising:	No
Yes	· ·	
Hygroscopic:	Irritant:	No
Yes		
Flammable:	Handling:	See caution, Section 2
No	-	
Other (specify): Contains material of human origin		
Toxicological properties		
Effects of inhalation: Not established, avoid inhalation		
Effects of ingestion: Not established, avoid ingestion		
Effects of skin absorption: Not established, avoid contact with skin		
Suggested First Aid		
Inhalation: Se	Seek medical advice	
Ingestion: Seek medical advice		
Contact with eyes: Wash with copious amounts of water. Seek		
medical advice		
Contact with skin: Wash thoroughly with water.		
Action on Spillage and Method of Disposal		

Spillage of ampoule contents should be taken up with absorbent material wetted with an appropriate disinfectant. Rinse area with an appropriate disinfectant followed by water.

Absorbent materials used to treat spillage should be treated as biological waste.

15. LIABILITY AND LOSS

In the event that this document is translated into another language, the English language version shall prevail in the event of any inconsistencies between the documents.

Unless expressly stated otherwise by NIBSC, NIBSC's Standard Terms and Conditions for the Supply of Materials (available at http://www.nibsc.org/About_Us/Terms_and_Conditions.aspx or upon request by the Recipient) ("Conditions") apply to the exclusion of all other terms and are hereby incorporated into this document by reference. The Recipient's attention is drawn in particular to the provisions of clause 11 of the Conditions.

Potters Bar, Hertfordshire, EN6 3QG. T +44 (0)1707 641000, nibsc.org WHO International Laboratory for Biological Standards, UK Official Medicines Control Laboratory