

WHO International Standard
Follicle Stimulating Hormone, (FSH, Urofollitropin), Urinary,
Human For Bioassay
NIBSC code: 92/512
Instructions for use
(Version 3.0, Dated 22/01/2008)

1. INTENDED USE

This consists of a batch of ampoules (coded 92/512) containing highly purified human urinary follicle stimulating hormone obtained from menopausal urine. It was established as the First International Standard for Follicle Stimulating Hormone, Urinary, Human for bioassay by the Expert Committee on Biological Standardization of the World Health Organisation at its forty sixth meeting in 1995 (1). It is intended for use in calibration of preparations of urinary FSH by bioassay.

2. CAUTION

This preparation is not for administration to humans or animals in the human food chain.

The preparation contains material of human origin, and either the final product or the source materials, from which it is derived, have been tested and found negative for HBsAg, anti-HIV and HCV RNA. As with all materials of biological origin, this preparation should be regarded as potentially hazardous to health. It should be used and discarded according to your own laboratory's safety procedures. Such safety procedures should include the wearing of protective gloves and avoiding the generation of aerosols. Care should be exercised in opening ampoules or vials, to avoid cuts.

3. UNITAGE

This preparation has been assigned a unitage of 121 International Units per ampoule.

4. CONTENTS

Country of origin of biological material: United Kingdom. Each ampoule contains the residue, after freeze drying, of 1ml of a solution which contained:

FSH approximately 10 µg
Human plasma albumin 2mg
Mannitol 10mg
Sodium chloride 15.4mM

5. STORAGE

Unopened ampoules should be stored at -20°C

Please note: because of the inherent stability of lyophilized material, NIBSC may ship these materials at ambient temperature.

6. DIRECTIONS FOR OPENING

DIN ampoules have an 'easy-open' coloured stress point, where the narrow ampoule stem joins the wider ampoule body. Various types of ampoule breaker are available commercially. To open the ampoule, tap the ampoule gently to collect material at the bottom (labelled) end and follow manufactures instructions provided with the ampoule breaker.

7. USE OF MATERIAL

For practical purposes each ampoule contains the same amount of the same materials. Dissolve all the contents in a known amount of buffer

solution. No attempt should be made to weigh portions of the freeze dried powder.

The material has not been sterilised and contains no bacteriostat.

8. STABILITY

It is the policy of WHO not to assign an expiry date to their international reference materials. They remain valid with the assigned potency and status until withdrawn or amended.

Reference materials are held at NIBSC within assured, temperature-controlled storage facilities. Reference Materials should be stored on receipt as indicated on the label. For information specific to a particular biological standard, contact standards@nibsc.ac.uk.

In addition, once reconstituted, diluted or aliquoted, users should determine the stability of the material according to their own method of preparation, storage and use.

NIBSC follows the policy of WHO with respect to its reference materials.

Users who have data supporting any deterioration in the characteristics of any reference preparation are encouraged to contact NIBSC.

9. REFERENCES

1. WHO Tech Rep. Series No.872, 1988

Further information can be obtained as follows;

http://www.nibsc.org/terms_and_conditions.aspx

10. FURTHER INFORMATION

This material: enquiries@nibsc.org
WHO Biological Standards:
http://www.who.int/biologicals/en/
JCTLM Higher order reference materials:
http://www.bipm.org/en/committees/jc/jctlm/
Derivation of International Units:
http://www.nibsc.org/standardisation/international_standards.aspx
Ordering standards from NIBSC:
http://www.nibsc.org/products/ordering.aspx
NIBSC Terms & Conditions:

11. CUSTOMER FEEDBACK

Customers are encouraged to provide feedback on the suitability or use of the material provided or other aspects of our service. Please send any comments to enquiries@nibsc.org

12. CITATION

In all publications, including data sheets, in which this material is referenced, it is important that the preparation's title, its status, the NIBSC code number, and the name and address of NIBSC are cited and cited correctly.





NIBSC Confidence in Biological Medicines

13. MATERIAL SAFETY SHEET

Classification in accordance with Directive 2000/54/EC, Regulation (EC) No 1272/2008: Not applicable or not classified

Physical and Chemical properties		
Friysical and Chemical properties		
Physical	Corrosive:	No
appearance: Freeze		
dried powder		
Stable:	Oxidising:	No
Yes		
Hygroscopic:	Irritant:	No
Yes	Literatilla au	Occupation Occident
Flammable:	Handling:	See caution, Section 2
No Other (anasity)		of house a spinia. Can page to the
Other (specify): Contains material of human origin. Can react with oxidising materials. Avoid contact with acids and alkalis.		
Oxidising materials. Avoid contact with acids and alkalis.		
Toxicological properties		
Effects of inhalation: Not established, avoid inhalation		
Effects of ingestion: Not established, avoid ingestion		
Effects of skin absorption: Not established, avoid contact with skin		
Suggested First Aid		
Inhalation: Seek medical advice		
Ingestion: Seek medical advice		
Contact with eyes: W	ash with copio	us amounts of water. Seek
medical advice		
Contact with skin: W	ash thoroughly with water.	
Action on Spillage and Method of Disposal		

14. LIABILITY AND LOSS

biological waste.

appropriate disinfectant followed by water.

In the event that this document is translated into another language, the English language version shall prevail in the event of any inconsistencies between the documents.

Spillage of ampoule contents should be taken up with absorbent material wetted with an appropriate disinfectant. Rinse area with an

Absorbent materials used to treat spillage should be treated as

Unless expressly stated otherwise by NIBSC, NIBSC's Standard Terms and Conditions for the Supply of Materials (available at http://www.nibsc.org/About_Us/Terms_and_Conditions.aspx or upon request by the Recipient) ("Conditions") apply to the exclusion of all other terms and are hereby incorporated into this document by reference. The Recipient's attention is drawn in particular to the provisions of clause 11 of the Conditions.

15. INFORMATION FOR CUSTOMS USE ONLY

Country of origin for customs purposes*: United Kingdom
* Defined as the country where the goods have been produced and/or
sufficiently processed to be classed as originating from the country of
supply, for example a change of state such as freeze-drying.

Net weight: 13mg
Toxicity Statement: Non-toxic

Veterinary certificate or other statement if applicable.

Attached: No

17. CERTIFICATE OF ANALYSIS

NIBSC does not provide a Certificate of Analysis for WHO Biological Reference Materials because they are internationally recognised primary reference materials fully described in the instructions for use. The reference materials are established according to the WHO Recommendations for the preparation, characterization and

National Institute for Biological Standards and Control,

Potters Bar, Hertfordshire, EN6 3QG. T +44 (0)1707 641000, nibsc.org WHO International Laboratory for Biological Standards, UK Official Medicines Control Laboratory

establishment of international and other biological reference standards http://www.who.int/bloodproducts/publications/TRS932Annex2_Inter_biol efstandardsrev2004.pdf (revised 2004). They are officially endorsed by the WHO Expert Committee on Biological Standardization (ECBS) based on the report of the international collaborative study which established their suitability for the intended use.

