

Non WHO Reference Material Rec. Staphylokinase (STAR) NIBSC Reagent NIBSC code: 94/718 Instructions for use (Version 6.0, Dated 04/01/2011)

This material is not for in vitro diagnostic use.

1. INTENDED USE

This preparation is intended for use as a standard in plasminogen activation assays designed to measure the potency of staphylokinase solutions

2. CAUTION

This preparation is not for administration to humans or animals in the human food chain.

The preparation contains material of human origin, and either the final product or the source materials, from which it is derived, have been tested and found negative for HBsAg, anti-HIV and HCV RNA. As with all materials of biological origin, this preparation should be regarded as potentially hazardous to health. It should be used and discarded according to your own laboratory's safety procedures. Such safety procedures should include the wearing of protective gloves and avoiding the generation of aerosols. Care should be exercised in opening ampoules or vials, to avoid cuts.

3. UNITAGE

Each ampoule contains 1.0 unit of recombinant Staphylokinase (STAR). The assigned content was valid at the time of manufacture-no information available on long term stability.

4. CONTENTS

Country of origin of biological material: United Kingdom.

The source material was supplied as a frozen solution of 25ml of STAR concentrate (1.05mg/ml) in a solution of 15mM potassium phosphate, 150mM NaCl, pH 7.4. This frozen solution was thawed and diluted with a 0.01M phosphate buffer (pH 7.4) containing 5mgs of human albumin per ml to a final concentration of 15 micrograms of STAR per ml. This batch of ampoules, coded 94/718, contains the residue after freeze drying of an accurately measured 1ml of a Staphylokinase solution (15 micrograms STAR in 10mM phosphate buffer, pH7.4 containing 5mg human albumin).

5. STORAGE

Unopened ampoules should be stored in the dark at or below –20°C. Please note: because of the inherent stability of lyophilized material, NIBSC may ship these materials at ambient temperature.

6. DIRECTIONS FOR OPENING

DIN ampoules have an 'easy-open' coloured stress point, where the narrow ampoule stem joins the wider ampoule body. Various types of ampoule breaker are available commercially. To open the ampoule, tap the ampoule gently to collect material at the bottom (labelled) end and follow manufactures instructions provided with the ampoule breaker.

7. USE OF MATERIAL

No attempt should be made to weigh out any portion of the freeze-dried material prior to reconstitution

The total contents of the ampoule should be reconstituted with 1.0ml distilled water and dissolved gently by swirling to avoid froth. The reconstituted reagent should be used as soon as possible after reconstitution

8. STABILITY

Reference materials are held at NIBSC within assured, temperaturecontrolled storage facilities. Reference Materials should be stored on receipt as indicated on the label.

NIBSC follows the policy of WHO with respect to its reference materials.

Following ampouling by procedures outlined by Campbell (1974) the stability was established by assaying the contents of ampoules stored at elevated temperatures (Jerne & Perry, 1956). It was observed that the reagent was stable for 18 months at 37°C using a modified clot lysis assay (Heath & Gaffney, 1990).

9. REFERENCES

Campbell, P.J., (1974) J. Biol. Stand. 2: 259-267

Heath, A. Gaffney, P.J, (1990) Thromb. Haemost. 64, 267-269

Jerne, N.K, Perry, W.L.M. (1956) Bull. WHO 14: 167-182.

10. ACKNOWLEDGEMENTS

We are grateful to the Hans Knoll Institute and the good offices of Professor D. Collen, for supplying the staphylokinase source material.

11. FURTHER INFORMATION

Further information can be obtained as follows; This material: enquiries@nibsc.org WHO Biological Standards: http://www.who.int/biologicals/en/ JCTLM Higher order reference materials: http://www.bipm.org/en/committees/jc/jctlm/ Derivation of International Units: http://www.nibsc.org/standardisation/international_standards.aspx Ordering standards from NIBSC: http://www.nibsc.org/products/ordering.aspx NIBSC Terms & Conditions: http://www.nibsc.org/terms_and_conditions.aspx

12. CUSTOMER FEEDBACK

Customers are encouraged to provide feedback on the suitability or use of the material provided or other aspects of our service. Please send any comments to enquiries@nibsc.org

13. CITATION

In all publications, including data sheets, in which this material is referenced, it is important that the preparation's title, its status, the NIBSC code number, and the name and address of NIBSC are cited and cited correctly.

14. MATERIAL SAFETY SHEET

Classification in accordance with Directive 2000/54/EC, Regulation (EC) No 1272/2008: Not applicable or not classified

| Physical and Chemical properties | | | |
|--|------------|----|--|
| Physical appearance: Freeze- dried powder | Corrosive: | No | |
| Stable: Yes | Oxidising: | No | |
| Hygroscopic: Yes | Irritant: | No | |

National Institute for Biological Standards and Control,

Potters Bar, Hertfordshire, EN6 3QG. T +44 (0)1707 641000, nibsc.org WHO International Laboratory for Biological Standards,

UK Official Medicines Control Laboratory

Medicines & Healthcare products Regulatory Agency



| Flammable: | Handling: | See caution, Section 2 | | | |
|--|--|--------------------------|--|--|--|
| No | 5 | | | | |
| Other (specify): Contains material of human origin | | | | | |
| Other (specify). Contains material of numari origin | | | | | |
| Toxicological properties | | | | | |
| | | | | | |
| Effects of inhalation: Not established, avoid inhalation | | | | | |
| Effects of ingestion: Not established, avoid ingestion | | | | | |
| Effects of skin absorption: Not established, avoid contact with skin | | | | | |
| | • · · · · | | | | |
| Suggested First Aid | | | | | |
| Inhalation: | Seek medical adv | ice | | | |
| Ingestion: Seek medical advice | | | | | |
| Contact with eyes: | Nash with copious | s amounts of water. Seek | | | |
| medical advice | • | | | | |
| Contact with skin: | Nash thoroughly | with water. | | | |
| | | | | | |
| Action on Spillage and Method of Disposal | | | | | |
| Spillage of ampoule contents should be taken up with absorbent | | | | | |
| | material wetted with an appropriate disinfectant. Rinse area with an | | | | |

material wetted with an appropriate disinfectant. Rinse area with an appropriate disinfectant followed by water. Absorbent materials used to treat spillage should be treated as

biological waste.

15. LIABILITY AND LOSS

In the event that this document is translated into another language, the English language version shall prevail in the event of any inconsistencies between the documents.

Unless expressly stated otherwise by NIBSC, NIBSC's Standard Terms and Conditions for the Supply of Materials (available at http://www.nibsc.org/About_Us/Terms_and_Conditions.aspx or upon request by the Recipient) ("Conditions") apply to the exclusion of all other terms and are hereby incorporated into this document by reference. The Recipient's attention is drawn in particular to the provisions of clause 11 of the Conditions.

16. INFORMATION FOR CUSTOMS USE ONLY

Country of origin for customs purposes*: United Kingdom * Defined as the country where the goods have been produced and/or sufficiently processed to be classed as originating from the country of supply, for example a change of state such as freeze-drying. Net weight: 10mg Toxicity Statement: Toxicity not assessed Veterinary certificate or other statement if applicable.

Attached: No

National Institute for Biological Standards and Control, Potters Bar, Hertfordshire, EN6 3QG. T +44 (0)1707 641000, nibsc.org WHO International Laboratory for Biological Standards, UK Official Medicines Control Laboratory