


NR-45833 Synfluenza (Synthetic Influenza) Clone Set, Recombinant in Escherichia coli, Plate 15 (Neuraminidase) (Clones)

[下载为PDF](#)

• 2次围观

产品图片



产品英文名称

[NR-45833_Synfluenza \(Synthetic Influenza\) Clone Set, Recombinant in Escherichia coli, Plate 15 \(Neuraminidase\)\(Clones\)](#)

产品别名

[NR-45833_Synfluenza \(Synthetic Influenza\) Clone Set, Recombinant in Escherichia coli, Plate 15 \(Neuraminidase\)\(Clones\)](#)

[NR-45833 Synfluenza \(Synthetic Influenza\) Clone Set, Recombinant in Escherichia coli, Plate 15 \(Neuraminidase\) \(Clones\)](#)

货号/SKU

NR-45833

货号/规格

EA

库存与交货期

3-8周

人民币价格

14000

人民币价格说明

本商品人民币2024年销售价格正在调整中，请等待更新完毕。

本商品的展示人民币价格已包含商品本身全部、VAT市担税13%、国际运输运费、国内物流运费、运输保险，以及冷试剂包装材料（例如液氮罐、泡沫箱、金属箱、蓝冰、湿冰、干冰、蓄冷剂、液氮等）、装卸费、报关资料费、人力支出等一切费用。

本商品的美元价、市场价、零售价、厂商指导价或该商品的曾经展示过的销售价等，并非商品原价，仅供参考。

试剂海关审批

A/B级风险物质只能直接使用者购买并持有实验室有效资质，其它询客服确认

国外采购

支持/部分限制一年内购买数量

厂牌

BEI Resources(ATCC)

品牌

[BEI Resources](#)

产品基础信息

生物安全等级建议分类：美国、1

产品描述信息

NR-45833?? Synfluenza (Synthetic Influenza) Clone Set, Recombinant in Escherichia coli, Plate 15 (Neuraminidase)(Clones)|Influenza A virus|Synfluenza (Synthetic Influenza) Clone Set, Recombinant in Escherichia coli, Plate 15 (Neuraminidase)-60°C or colder|Pathogen Functional Genomics Resource Center at the J. Craig Venter Institute|Acknowledgment for publications should read "The following reagent was obtained through BEI Resources, NIAID, NIH: Synfluenza (Synthetic Influenza) Clone Set, Recombinant in *Escherichia coli*, Plate 15 (Neuraminidase), NR-45833."|Quantity limit per order for this item is 1. This item can be ordered twice a year. Orders over this limit will be sent to NIAID for approval before shipment.

This item is currently in our production queue. Please allow ample time for distribution lots to be made available.

The Synfluenza clone set is part of a National Institute of Allergy and Infectious Diseases (NIAID) initiative to create 1000 influenza gene segment clones from 12 host subtypes that span the protein sequence diversity of influenza viruses between 2005 and 2010. Each clone is designed from GenBank sequences with consensus untranslated regions. The purpose of the project is to develop the ability to create and stockpile synthetic DNA encoding influenza gene segments. These segments can then be used to generate virus seed stocks and a library of clones for vaccine, diagnostic and basic research.

The NIAID Genome Sequencing Center at the J. Craig Venter Institute constructed synthetic influenza neuraminidase (NA) and hemagglutinin (HA) genes using automated DNA synthesis and assembly. There are nine synthetic NA influenza clone plates (BEI numbers NR-45827 through NR-45833, NR-45090 and NR-45091) and six synthetic HA influenza clone plates (BEI numbers NR-45092 through NR-45097) in the set.

Each synthetic NA gene from NR-45833 was manufactured from five individually-designed, double-stranded DNA construct cassettes produced by assembly of eight chemically-synthesized oligonucleotides using the Gibson Assembly[®] process. The five cassettes were combined into the pSMART[™]-LCKan vector (Lucigen[™]) to establish gene segment clones in One Shot[™] TOP10 competent (Invitrogen[™]) *Escherichia coli* (*E. coli*) cells.

Detailed information provided by JCVI for the complete Synfluenza Clone Set, including insert sequence, is available in the attached Synfluenza Master Clone List.

Each well of the 96-well plate contains approximately 200 μ L of *E. coli* culture in Yeast Extract Tryptone media containing 25 μ g/mL kanamycin supplemented with 10% glycerol.

Note: Production in the 96-well format has increased risk of cross-contamination between adjacent wells. Individual clones should be purified (e.g. single colony isolation and purification using good microbiological practices) and sequence-verified prior to use.

主要内容

此项目的每个订单数量限制为1.此商品每年可订购两次.此限制的订单将被发送到发货前的批准.此项目目前在我们的生产队列中. nbsp;nbsp;nbsp;请允许有充分的时间进行分销批次. Synfluenza Clone Set是国家过敏和传染病(NIAID)倡议的一部分,从12月开始创造1000个流感基因克隆培养2005年至2010年间流感病毒蛋白序列多样性的宿主亚型.每个克隆都是由Genbank序列的设计,不合适的区域.该项目的目的是制定制造编码流感基因段的合成DNA的能力.然后可以使用这些段来产生病毒种子股和克隆文库,用于疫苗,诊断和基础研究. J.Craig Venter Institute建造合成流感神经氨酸酶(NA)的NIAID基因组测序中心血凝素(HA)基因使用自动DNA合成和组装.含有九种合成Na流感克隆平板(BEI编号NR-45827至NR-45833, NR-45090和NR-45091)和六种合成的HA流感克隆平板(BEI编号NR-45092至NR-45097).来自NR-45833的每个合成NA基因由五个单独设计的双链DNA构建盒制成,通过组装八种化学合成的寡核苷酸使用Gibson组装?方法制备.将五个盒式磁带组合成PSMART?-Lkan载体(Lucigen?),以在一个镜头中建立基因克隆? top10 possent (Invitrogen?)大肠杆菌(大肠杆菌)细胞. JCVI为完整的Synfluenza Clone集,包括插入序列,附加的Synfluenza Master Clone列表提供了详细信息. 96孔板的每个孔含有约200 μ l的 e. Coli 在酵母提取物中的培养蛋白培养基,其含有25 μ g/ml卡那霉素,其补充有10%甘油. 注意: 96孔格式的生产增加了风险增加相邻孔之间的交叉污染.应纯化各个克隆(例如,使用良好的微生物学实践的分离和纯化)和在使用前进行序列验证.

厂牌介绍

BEI Resources 由美国国家过敏和传染病研究所(NIAID)成立,旨在为研究A、B和C类优先病原体、**新兴传染病**病原体、非病原微生物和其他相关微生物材料提供试剂、工具和材料到研究界. BEI Resources 获取、验证和生产科学家进行基础研究和开发改进的诊断测试、疫苗和疗法所需的试剂. 通过将这些功能集中在 BEI Resources 中,可以监控科学界对这些材料的访问和使用,并确保试剂的质量控制.

除了为传染病界提供材料外, BEI Resources 还鼓励和支持研究人员和机构的材料存放. **使用 BEI Resources 存放材料**对研究人员和研究社区有许多优势,包括安全存储、社区访问和分发;同时保护存款人的知识产权. 只要有需要, BEI 资源库将作为研究人员的资源进行维护. 您在 BEI Resources 的存款是一项有助于未来研究的长期投资.

BEI Resources 自 2003 年起由**美国典型培养物保藏中心(ATCC) 根据合同管理. 2016 年 5 月, ATCC**获得了一份为期七年的继续管理 BEI Resources 的合同,合同范围已扩大到更全面的研究目录材料,包括由其他政府支持的研究项目存放的材料,将提供给生物防御和新兴传染病科学界. 真菌、寄生虫、载体和其他相关材料已添加到现有的细菌、病毒和毒素试剂中,涵盖 NIAID A、B 和 C 类优先病原体和 NIAID 指定的新发传染病病原体和生物.

品牌标识

beii RESOURCES

SUPPORTING INFECTIOUS DISEASE RESEARCH

产品关键词

[BEI Resources 生物材料代理进口报关公司](#)

[如何购买可以直接从 BEI Resources 购买菌株吗](#)

[BEI Resources 生物材料库中国官网是?](#)

[BEI Resources 生物材料库中国代理](#)

[BEI Resources 怎么查询菌株抗体细胞](#)

[NR-45833 Synfluenza \(Synthetic Influenza\) Clone Set](#)

[Recombinant in Escherichia coli](#)

[Plate 15 \(Neuraminidase\) \(Clones\), NR-1 Vaccinia virus](#)

[Modified Vaccinia Ankara \(MVA\)\(Viruses\)](#)

[Modified Vaccinia Ankara \(MVA\), 生物试剂报关 BEI Resources](#)

[NIAID](#)

[NIH](#)

[NIAID 全称是美国国家过敏和传染病研究所](#)

[一键获取大包装优惠价](#)

无

选择您的报价场景

- 【我们直接使用】需要优惠价、大包装规格、货期 -- -- --> [报价默认含增值税13%发票; 尽量提供货号、规格、需求数量]
- 【需要技术文档】产品说明书、COA、MSDS、手册 -- -- --> [默认提供说明书或者COA, 特别技术指标要求请下面填入详细描述]
- 【我帮客户找货】需要优惠价、大包装规格、货期 -- -- --> [报价默认含增值税13%发票]
- 【推荐替代产品】需要优惠价、大包装规格、货期 -- -- --> [提供替代产品的价格, 默认含增值税13%发票]
- 【我能原厂直采】请只提供代理进口清关服务的报价 -- -- --> [适合只需要进口许可证代办服务、清关服务的专业级买家, 独立服务]
- 【其它报价场景】

请输入您的情况与报价要求

报价要求详细描述

[如有请填写;若不留空即可] 按10KG、25L大量采购的时候, 是否可

贵单位贵姓

接受报价的E-mail

[马上发我报价](#)

手机扫描二维码阅读本页



可能感兴趣的内容

[Tris-HCl缓冲液\(1.5mol/L, pH8.8\) / 1.5M Tris-HCl缓冲液, pH8.8 \(BZ201\) 1000ml](#)

2021-12-13

[NR-46375小葡萄球菌_VCU012\(NRS847\)\(细菌\)](#)

2022-04-01

[NR-48263金黄色葡萄球菌亚种.金黄色葡萄球菌_USA300JE2.转座子突变体SAUSA300_1203\(NE1721\)\(突变细菌\)](#)

2022-04-01

[NR-2497嗜热脂肪地芽孢杆菌_NCA1518\(细菌\)](#)

2022-03-31

[甲状腺激素_人型_合成\[1-34\]片段82/508](#)

2024-05-19

[北京大学彩虹实验室招聘博士后](#)

2022-01-10

[水凝胶掩膜板精密切割pet电极片加工PDMS薄膜麦拉片异形激光裁切](#)

2021-12-02

[NR-48042金黄色葡萄球菌亚种.金黄色葡萄球菌_USA300JE2.转座子突变体SAUSA300_0043\(NE1500\)\(突变细菌\)](#)

2022-04-01

[NR-45989金黄色葡萄球菌_A970627\(细菌\)](#)

2022-04-01

[阿尔瓦雷斯龙快速小型化和适应性辐射研究获进展](#)

2020-08-04

[NR-46360金黄色葡萄球菌_OR-589\(NRS834\)\(细菌\)](#)

2022-04-01

[尹玉新课题组发现单链DNA结合蛋白RPA1蛋白全新的生理功能](#)

2022-06-17

[NR-48939免疫球蛋白,3476\(真菌\)](#)

2022-04-01

[基于改进YOLOv5的船舶多尺度SAR图像检测算法](#)

2024-02-26

[人造模拟合成眼泪、合成泪液 \(BZ182\) 1000ml](#)

2021-12-13

[TRIS乙酸 \(TAE缓冲\) 运行缓冲液 \(50xTAE\) \(BZ208\) 200ml](#)

2021-12-13

[NR-46652金黄色葡萄球菌亚种,金黄色葡萄球菌,USA300JE2,转座子突变体SAUSA300_1668\(NE109\)\(突变细菌\)](#)

2022-04-01

[DIN EN 1616:1999人工模拟合成尿液 \(BZ102\) 500ml](#)

2021-12-13

[香蕉种质资源鉴定获得了新进展](#)

2022-06-17

[破伤风抗毒素,马,用于生物测定60/013](#)

2024-05-19