

[D4265-14人工模拟合成皮脂液体-无防腐剂（BZ116）500ml](#)

[下载为PDF](#)

- 17 次围观

产品图片



产品英文名称

[D4265-14 Artificial Sebum - Not Stabilized \(BZ116\)](#)

产品别名

[人工模拟体液](#)

货号/SKU

Chemazone306

库存与交货期

1-2周

人民币价格

14493

人民币价格说明

本商品人民币2024年销售价格正在调整中，请等待更新完毕。

本商品的展示的人民币价格已包含商品本身金额、VAT增值税13%、国际运输运费、国内物流运费、运输保险、以及冷链包装材料（例如液氮罐、泡沫箱、金属桶、蓝冰、湿冰、干冰、蓄冷剂、液氮等）、装卸费、相关资料费、人力支出等一切费用。

试剂海关审批

支持

国外采购

支持

厂牌

Chemazone Inc

品牌

[Biochemazone™](#)

产品基础信息

1000mL, 100mL, 200mL, 500mL

产品描述信息

Artificial Sebum

Sebum is an oily secretion produced by sebaceous glands, which spreads over the hair and skin for waterproofing purposes. Artificial Sebum is used for the study of home laundry products and conditions to remove stains from fabric. Artificial [eccrine](#) sebum for CORROSION RESISTANCE and Alloy study applications. Buy Premium quality Simulated [Sweat](#) Sebum solution at Low price.

Artificial Sebum or Simulated Sebum (Custom pH and stabilized versions also available)

This reagent is only used in the field of scientific research and for in vitro use only.

Pack Size: 100 mL, 200 mL, 500 mL, 1000 mL

Validity: 6 months for Not Stabilized and up to 12 months for Stabilized version. At -20 Deg C, in the Frozen stage, shelf life is up to 1 year for non-stabilized versions.

Transportation and storage: transportation at room temperature, storage at 2-8 °C.

Artificial Sebum Features & Note:

1. Simulated Sebum is a sterile solution and can be used directly. Our Artificial Sebum provides the reliability, reproducibility, and convenience needed for testing.
2. Simulated Sebum, store at 4 degrees, transported at room temperature.
3. This Sebum emulsion is limited to the scientific research of professionals, and cannot be used for clinical diagnosis or treatment, and cannot be used for food or medicine, and cannot be stored in ordinary houses.
4. For your safety and health, please wear a lab coat and disposable gloves.

Biochemazone™ provides an [artificial sebum](#) formulation according to ASTM designation D4265-14 or D4265-98. This formulation can be used according to the AATCC Standard Test Method 130-2010.

[Artificial Sebum](#) consists of amino acids, minerals, and various metabolites such as Oleic acid as an emulsifying agent which simulates the composition, properties, and pH of the real Human Sebum. Buy [eccrine sweat](#) at Low price. We add preservatives, stabilizers, antifungal, and antibacterial reagents to enhance the shelf life. Please request it while making the order.

[Biochemazone™](#) offers a diverse selection of [Artificial fluids](#) for the testing & Diagnostic of various diseases in Vitro laboratory and Research. Check out our various state biofluid inventory to see what's ready to ship today.

For a custom request, please send us an email at sales@coreab.cn

产品安全信息

以下中文仅供参考，如专业术语有误，请以英文为准！###人造皮脂 皮脂是通过皮脂腺产生的油性分泌物，其在头发和皮肤上延伸以防防水目的。人造皮脂用于研究家庭洗衣产品和条件，以从织物中去除污渍。用于耐腐蚀性和合金研究应用的人工生态生态带。购买优质质量模拟汗水SEBUM解决方案，价格低廉。人造皮脂或模拟皮脂（也可用的定制pH和稳定版本）该试剂仅用于科学研究领域，仅用于体外使用。包装尺寸：100毫升，200毫升，500mL，1000毫升 有效性：稳定版不稳定，6个月，最多可达12个月。在-20°C下，在冷冻阶段，保质期长达1年的非稳定版本。运输和储存：室温运输，2-8°C储存。人工皮脂特征笔记：1.模拟皮脂是一种无菌溶液，可直接使用。我们的人工表皮可以提供测试所需的可靠性，再现性和便利性。2.模拟皮脂，储存在4度，在室温下运输。3.这种皮脂乳液仅限于对专业人士的科学研究，不能用于临床诊断或治疗，不能用于食物或药物，不能储存在普通房屋中。4.为了您的安全和健康，请戴上实验室外套和一次性手套。Biochemazone™根据ASTM指定D4265-14或D4265-98提供人造皮脂制剂。该配方可根据AATCC标准测试方法130-2010使用。人造皮脂由氨基酸，矿物质和各种代谢物如油酸，作为乳化剂，其模拟了真实人体皮脂的组成，性质和pH。以低廉的价格购买经济型汗水。我们添加防腐剂，稳定剂，抗真菌和抗菌试剂，以增强保质期。请在下单时请求。Biochemazone™为测试&提供提供各种人造液。各种疾病在体外实验室和研究的诊断。

查看我们各种国家生物流体库存，看看今天准备发货。有关自定义请求，[请向我们发送电子邮件至sales@coreab.cn](mailto:sales@coreab.cn)
[人造皮脂](#) [人造皮脂装备清单](#) [人造皮脂](#) [其他组合物和客户](#)；S指定的配置可以根据要求提供。稳定剂，pH值，包装尺寸和模拟皮脂流体的含量是可定制的。联系我们进行定制：Sales@coreab.cn 运输和储存：室温下的运输，4°C储存。对于长期储存在-20°C下冻结它 [/ vc_column_text]

主要内容

ARTIFICIAL SEBUM

ARTIFICIAL SEBUM LIST OF INGREDIENTS

ARTIFICIAL SEBUM

Other compositions and customer's specified configurations can be delivered upon request. Stabilizer, pH Value, Pack Size, and Contents of Simulated Sebum Fluid are customizable.

Contact us for Customization: sales@coreab.cn

Transportation and storage: Transportation at room temperature, storage at 4 °C. For Longterm Storage Freeze it at -20 °C

[/vc_column_text]

厂牌介绍

关于 Biochemazone™

人工模拟生物体液试剂和仪器的先驱

Biochemazone™是一家不断发展的生物科技公司，提供高品质、纯化的人工生物化学品、试剂和解决方案，如模拟唾液、模拟尿液、模拟汗液等，并开发学术研究机构、大学、医院、工业、政府和私人实验室。

我们的生化试剂经过全面测试，可满足对卓越品质、不妥协生物的研发需求，以具有竞争力的价格准确模拟真正的生物解决方案活动。

Biochemazone™在科学界建立了良好的声誉，其声誉建立在创新产品、客户服务和对细节的关注上。我们的质量声誉是在市场上取得有利地位以及与制造顶级合成生物产品（如人造唾液、人造尿液、人造皮脂、人造汗液、人造汗液、人造血清、人工血液、人工精液、人工耳垢等等.....

欢迎消费品测试实验室和制造商联系 **Biochemazone™**或本地代理商以了解人工模拟生物流体液体和生化试剂的需求。

品牌标识



产品关键词

[Biochemazone人造体液](#)

[Biochemazone™人造体液目录价格](#)

[Biochemazone人造汗液](#)

[人工汗液厂家](#)

[人造体液](#)

[模拟人体体液biochemazone](#)

[体液包括哪些](#)

[人体体液是什么](#)

[什么是人体体液](#)

[人造模拟唾液是体液吗](#)

[人体体液占体重多少](#)

[人造唾液](#)

[人造肺泡灌洗液](#)

[人工痰液](#)

[人造粪便](#)

[人工尿液](#)

一键获取大包装优惠报价

- 无 -

选择您的报价场景

- 【我们直接使用】需要优惠报价、大包装规格、货期 -- ---->[报价默认含增值税13%发票；尽量提供货号、规格、需求数量]
- 【需要技术文档】产品说明书、COA、MSDS、手册 -- ---->[默认提供说明书或者COA，特别技术指标要求请下面填入详细描述]
- 【我帮客户找货】需要优惠报价、大包装规格、货期 -- ---->[报价默认含增值税13%发票]
- 【推荐替代产品】需要优惠报价、大包装规格、货期 -- ---->[提供替代产品的价格，默认含增值税13%发票]
- 【我能原厂直采】请只提供代理进口清关服务的报价 -- ---->[适合只需要进口许可证代办服务、清关服务的专业级买家，独立服务]
- 【其它报价场景】

请输入您的情况与报价要求

报价要求详细描述

【如有请填写;若无留空即可】按10KG、25L大量采购的时候，是否可

贵单位贵姓

接受报价的E-mail

手机扫描二维码阅读本页



可能感兴趣的内容

[HM-34柠檬酸杆菌属,30_2\(细菌\)](#)

2022-04-01

[道康宁DC184 SYLGARD 184灌封胶PDMS 110g](#)

2021-12-02

[NR-13472结核分枝杆菌,CDC1551TransposonMutant34\(MT1862,Rv1814\)\(突变细菌\)](#)

2022-04-01

[NR-22091_肠沙门氏菌,SerovarNewport,分离物S8949\(细菌\)](#)

2022-04-01

[专家提出噬菌体防御系统基因的影响应引起重视](#)

2022-06-17

[MRA-742chabaudichabaudi疟原虫,AS\(3CQ\)\(寄生原生动\)](#)

2022-04-01

[NR-28911来自成年雌性曼氏血吸虫的基因组DNA,菌株NMRI\(核酸\)](#)

2022-04-01

[NR-3228多克隆抗仙台病毒,\(抗血清,豚鼠\)\(多克隆抗血清\)](#)

2022-03-31

[NR-46819金黄色葡萄球菌亚种.金黄色葡萄球菌,USA300JE2,转座子突变体SAUSA300_0293\(NE276\)\(突变体细菌\)](#)

2022-04-01

[NR-42566_晚期免疫性乙型脑炎病毒抗血清\(多克隆抗血清\)](#)

2022-04-01

[呼肠孤病毒\(REO \) FITC共轭物1毫升](#)

2019-05-08

[抗CD63 \[ME491\]抗体](#)

2021-12-21

[NR-50234_寨卡病毒,PLCa1_ZV?\(人类/2013/泰国\)\(病毒\)](#)

2022-04-01

[研究揭示隐藻的光适应与捕光调节机制](#)

2024-08-19

[NR-10133金黄色葡萄球菌,TCH1516\(细菌\)](#)

2022-04-01

[NR-3976多克隆抗黄斑支原体,PG-15,\(抗血清,驴\)\(多克隆抗血清\)](#)

2022-03-31

[NR-48130金黄色葡萄球菌亚种.金黄色葡萄球菌,USA300JE2,转座子突变体SAUSA300_1783\(NE1588\)\(突变细菌\)](#)
2022-04-01

[NR-29036_甲型流感病毒,A/PuertoRico/8-CV9/1934\(H1N1\)\(病毒\)](#)
2022-04-01

[NR-38044型登革热病毒,D85-019\(病毒\)](#)
2022-03-31

[NR-47063金黄色葡萄球菌亚种金黄色葡萄球菌、JE2、转座子突变体NE520\(SAUSA300_1452\)\(突变细菌\)](#)
2022-04-01