


# 重组抗体,抗VSV-G [8G5F11],人,IGMK

[下载为PDF](#)

- 5 次围观

产品图片



产品英文名称

[Recombinant Anti-VSV-G \[8G5F11\], Human, IgMk](#)

产品别名

[Kerafast独特的生物试剂](#)

货号/SKU

Kf-Ab01401-15.0

货号/规格

50ug

库存与交货期

4-5 周

人民币价格

13300

人民币价格说明

本产品人民币2024年销售价格正在调整中，请等待更新完毕。

本商品的展示的人民币价格已包含商品本身金额、VAT增值税13%、国际运输运费、国内物流运费、运输保险、以及冷链包装材料（例如液氮罐、泡沫箱、金属桶、蓝冰、湿冰、干冰、蓄冷剂、液氮等）、装卸费、相关资料费、人力支出等一切费用。

本商品的美元价、市场价、零售价、厂商指导价或该商品的曾经展示过的销售价等，并非商品原价，仅供参考。

试剂海关审批

使用人负责A/B风险申请资质

国外采购

支持/部分需签MTA

厂牌

Kerafast, Inc.

品牌

[Kerafast®](#)

产品基础信息

From the laboratory of Douglas S. Lyles, PhD, Wake Forest School of Medicine.

产品描述信息

Product Type:

Antibody

|                      |   |
|----------------------|---|
| Name:                | Anti-VSV-G [8G5F11]   |
| Antigen:             | VSV-G   |
| Host:                | Mouse   |
| Isotype:             | <b>EB0010:</b> IgG2a kappa<br><b>Recombinant versions:</b> see product name                             |
| Clonality:           | Monoclonal  |
| Clone Name:          | 8G5FL11 (I1)  |
| Specificity:         | VSV-Ind glycoprotein (G) protein  |
| Reactivity:          | Human   |
| Immunogen:           | VSV infection   |
| Format:              | Liquid  |
| Purity:              | Protein G purified  |
| Buffer:              | <b>EB0010:</b> PBS, 0.05% (w/v) Sodium Azide<br><b>Recombinant versions:</b> PBS with 0.02% Proclin 301 |
| Tested Applications: | Western blot (1:1000), flow cytometry, immunofluorescence   |
| Concentration:       | 1mg/mL  |
| Amount:              | 100uL   |
| Storage:             | -20C (avoid repeated freeze / thaw cycles)  |
| Shipped:             | Cold packs  |

#### 产品信息

Gregory DA, Olinger GY, Lucas TM, Johnson MC. Diverse viral glycoproteins as well as CD4 co-package into the same human immunodeficiency virus (HIV-1) particles. *Retrovirology*. 2014 Apr 3;11(1):28. Kubo S, Mitani K. A New Hybrid System Capable of Efficient Lentiviral Vector Production and Stable Gene Transfer Mediated by a Single Helper-Dependent Adenoviral Vector. *J Virol* 77: 2964-2971, 2003. Kobayashi M, Iida A, Ueda Y, and Hasegawa M. Pseudotyped Lentivirus Vectors Derived from Simian Immunodeficiency Virus SIVagm with Envelope Glycoproteins from Paramyxovirus. *J Virol* 77: 2607-2614, 2003. Burns JC, Friedmann T, Driever W, Burrascano M, Yee JK. Vesicular stomatitis virus G glycoprotein pseudotyped retroviral vectors: Concentration to very high titer and efficient gene transfer into mammalian and nonmammalian cells. *PNAS* 90: 8033-8037, 1993. Vandepol SB, Lefrancois L, Holland JJ. Sequences of the Major Antibody Binding Epitopes of the Indiana Serotype of Vesicular Stomatitis Virus. *Virology* 148: 312-325 (1986). Lefrancois L, Lyles DS. The interaction of antibody with the major surface glycoprotein of vesicular stomatitis virus. I. Analysis of neutralizing epitopes with monoclonal antibodies. *Virology* 121: 157-167, 1982. Boso G, Somia NV. Characterization of resistance to rhabdovirus and retrovirus infection in a human myeloid cell line. *PLoS One*. 2015 Mar 26;10(3):e0121455. doi: 10.1371/journal.pone.0121455. View Article Zhang J, Reiling C, Reinecke JB, Prislani I, Marky LA, Sorgen PL, Naslavsky N, Caplan S. Rabankyrin-5 interacts with EHD1 and Vps26 to regulate endocytic trafficking and retromer function. *Traffic*. 2012 May;13(5):745-57. View Article McNulty S, Flint M, Nichol ST, Spiropoulou CF. Host mTORC1 signaling regulates andes virus replication. *J Virol*. 2013 Jan;87(2):912-22. View Article Guzzo C, Fox J, Lin Y, Miao H, Cimbri R, Volkman BF, Fauci AS, Lusso P. The CD8-derived chemokine XCL1/lymphotactin is a conformation-dependent, broad-spectrum inhibitor of HIV-1. *PLoS Pathog*. 2013;9(12):e1003852. View Article Chen D, Gibson ES, Kennedy MJ. A light-triggered protein secretion system. *J Cell Biol*. 2013 May 13;201(4):631-40. View Article Xing M, Peterman MC, Davis RL, Oegema K, Shiao AK, Field SJ. GOLPH3 drives cell migration by promoting Golgi reorientation and directional trafficking to the leading edge. *Mol Biol Cell*. 2016 Dec 1;27(24):3828-3840. PubMed PMID: 27708138; PubMed Central PMCID: PMC5170606. View Article Yonemura Y, Li X, Müller K, Krämer A, Atigbire P, Mentrup T, Feuerhake T, Kroll T, Shomron O, Nohl R, Arndt HD, Hoischen C, Hemmerich P, Hirschberg K, Kaether C. Inhibition of cargo export at ER exit sites and the trans-Golgi network by the secretion inhibitor FLI-06. *J Cell Sci*. 2016 Oct 15;129(20):3868-3877. PubMed PMID: 27587840. View Article Yuan WC, Lee YR, Lin SY, Chang LY, Tan YP, Hung CC, Kuo JC, Liu CH, Lin MY, Xu M, Chen ZJ, Chen RH. K33-Linked Polyubiquitination of Coronin 7 by Cul3-KLHL20 Ubiquitin E3 Ligase Regulates Protein Trafficking. *Mol Cell*. 2014 May 22;54(4):586-600. doi: 10.1016/j.molcel.2014.03.035. PubMed PMID: 24768539. View Article McNulty S, Flint M, Nichol ST, Spiropoulou CF. Host mTORC1 signaling regulates andes virus replication. *J Virol*. 2013 Jan;87(2):912-22. doi: 10.1128/JVI.02415-12. PubMed PMID: 23135723; PubMed Central PMCID: PMC3554081. View Article Ferlin A, Raux H, Baquero E, Lepault J, Gaudin Y. Characterization of pH-sensitive molecular switches that trigger the structural transition of vesicular stomatitis virus glycoprotein from the postfusion state toward the prefusion state. *J Virol*. 2014 Nov;88(22):13396-409. doi: 10.1128/JVI.01962-14. PubMed PMID: 25210175; PubMed Central PMCID: PMC4249061. View Article Rao M, Song W, Jiang A, Shyr Y, Lev S, Greenstein D, Brantley-Sieders D, Chen J. VAMP-associated protein B (VAPB) promotes breast tumor growth by modulation of Akt activity. *PLoS One*. 2012;7(10):e46281. doi: 10.1371/journal.pone.0046281. PubMed PMID: 23049696; PubMed Central PMCID: PMC3462209. (supplemental information) View Article Swick A, Baltes A, Yin J. Visualizing infection spread: dual-color fluorescent reporting of virus-host interactions. *Biotechnol Bioeng*. 2014 Jun;111(6):1200-9. doi: 10.1002/bit.25170. PubMed PMID: 24338628; PubMed Central PMCID: PMC4004699. View Article Rogalin

HB, Heldwein EE. Characterization of Vesicular Stomatitis Virus Pseudotypes Bearing Essential Entry Glycoproteins gB, gD, gH, and gL of Herpes Simplex Virus 1. *J Virol*. 2016 Oct 28;90(22) View Article

Baquero E, Albertini AA, Raux H, Abou-Hamdan A, Boeri-Erba E, Ouldali M, Buonocore L, Rose JK, Lepault J, Bressanelli S, Gaudin Y. Structural intermediates in the fusion-associated transition of vesiculovirus glycoprotein. *EMBO J*. 2017 Mar 1;36(5):679-692. doi: 10.15252/embj.201694565. PubMed PMID: 28188244; PubMed Central PMCID: PMC5331758. View Article

Bowen AB, Bourke AM, Hiester BG, Hanus C, Kennedy MJ. Golgi-independent secretory trafficking through recycling endosomes in neuronal dendrites and spines. *Elife*. 2017 Sep 6;6. pii: e27362. View Article

Sakata M, Tani H, Anraku M, Kataoka M, Nagata N, Seki F, Tahara M, Otsuki N, Okamoto K, Takeda M, Mori Y. Analysis of VSV pseudotype virus infection mediated by rubella virus envelope proteins. *Sci Rep*. 2017 Sep 14;7(1):11607. View Article

Anandi L, Chakravarty V, Ashiq KA, Bodakuntla S, Lahiri M. DNA-PK plays a central role in transformation of breast epithelial cells following alkylating damage. *J Cell Sci*. 2017 Sep 18. pii: jcs.203034. View Article

Milani M, Annoni A, Bartolaccini S, Biffi M, Russo F, Di Tomaso T, Raimondi A, Lengler J, Holmes MC, Scheiflinger F, Lombardo A, Cantore A, Naldini L. Genome editing for scalable production of alloantigen-free lentiviral vectors for in vivo gene therapy. *EMBO Mol Med*. 2017 Nov;9(11):1558-1573. View Article

Tezuka K, Okuma K, Kuramitsu M, Matsuoka S, Tanaka R, Tanaka Y, Hamaguchi I. Control of HTLV-1 Infection by Eliminating Envelope Protein-Positive Cells with Recombinant Vesicular Stomatitis Viruses Encoding HTLV-1 Primary Receptor. *J Virol*. 2017 Dec 6. pii: JVI.01885-17. View Article

Nikolic J, Belot L, Raux H, Legrand P, Gaudin Y, Albertini A. Structural basis for the recognition of LDL-receptor family members by VSV glycoprotein. *Nat Commun*. 2018 Mar 12;9(1):1029. View Article

Felt SA, Droby GN, Grzelishvili VZ. Ruxolitinib and Polycation Combination Treatment Overcomes Multiple Mechanisms of Resistance of Pancreatic Cancer Cells to Oncolytic Vesicular Stomatitis Virus. *J Virol*. 2017 Jul 27;91(16). pii: e00461-17. View Article

Ci Y, Yang Y, Xu C, Shi L. Vesicular stomatitis virus G protein transmembrane region is crucial for the hemi-fusion to full fusion transition. *Sci Rep*. 2018 Jul 13;8(1):10669. View Article

Tijani M, Munis AM, Perry C, Sanber K, Ferrareso M, Mukhopadhyay T, Themis M, Nisoli I, Mattiuzzo G, Collins MK, Takeuchi Y. Lentivector Producer Cell Lines with Stably Expressed Vesiculovirus Envelopes. *Mol Ther Methods Clin Dev*. 2018 Aug 7;10:303-312. View Article

Johnson B, VanBlargan LA, Xu W, White JP, Shan C, Shi PY, Zhang R, Adhikari J, Gross ML, Leung DW, Diamond MS, Amarasinghe GK. Human IFIT3 Modulates IFIT1 RNA Binding Specificity and Protein Stability. *Immunity*. 2018 Mar 20;48(3):487-499.e5. View Article

Rahajeng J, Kuna RS, Makowski SL, Tran TTT, Buschman MD, Li S, Cheng N, Ng MM, Field SJ. Efficient Golgi Forward Trafficking Requires GOLPH3-Driven, PI4P-Dependent Membrane Curvature. *Dev Cell*. 2019 Jun 7. pii: S1534-5807(19)30451-4. View Article

Munis AM, Tijani M, Hassall M, Mattiuzzo G, Collins MK, Takeuchi Y. Characterization of Antibody Interactions with the G Protein of Vesicular Stomatitis Virus Indiana Strain and Other Vesiculovirus G Proteins. *J Virol*. 2018 Nov 12;92(23). pii: e00900-18. View Article

Moskovskich A, Goldmann U, Kartnig F, Lindinger S, Konecka J, Fiume G, Girardi E, Superti-Furga G. The transporters SLC35A1 and SLC30A1 play opposite roles in cell survival upon VSV virus infection. *Sci Rep*. 2019 Jul 18;9(1):10471. View Article

Gale TV, Horton TM, Hoffmann AR, Branco LM, Garry RF. Host Proteins Identified in Extracellular Viral Particles as Targets for Broad-Spectrum Antiviral Inhibitors. *J Proteome Res*. 2019 Jan 4;18(1):7-17. View Article

Moskovskich A, Goldmann U, Kartnig F, et al. The transporters SLC35A1 and SLC30A1 play opposite roles in cell survival upon VSV virus infection. *Sci Rep*. 2019;9(1):10471. Published 2019 Jul 18. View article

Rosor AM, Sleigh JN, Groves M, et al. Loss of BICD2 in muscle drives motor neuron loss in a developmental form of spinal muscular atrophy. *Acta Neuropathol Commun*. 2020;8(1):34. Published 2020 Mar 17. View article

Lucken-Ardjomande Häsler S, Vallis Y, Pasche M, McMahon HT. GRAF2, WDR44, and MICAL1 mediate Rab8/10/11-dependent export of E-cadherin, MMP14, and CFTR  $\Delta F508$ . *J Cell Biol*. 2020;219(5):e201811014. View article

Van Bergen NJ, Ahmed SM, Collins F, et al. Mutations in the exocyst component EXOC2 cause severe defects in human brain development. *J Exp Med*. 2020;217(10):e20192040. View article

Johnson MC, Lyddon TD, Suarez R, Salcedo B, LePique M, Graham M, Ricana C, Robinson C, Ritter DG. Optimized pseudotyping conditions for the SARS COV-2 Spike glycoprotein. *J Virol*. 2020 Aug 11;JVI.01062-20. View article

Azad T, Singaravelu R, Crupi MJF, Jamieson T, Dave J, Brown EEF, Rezaei R, Taha Z, Boulton S, Martin NT, Surendran A, Poutou J, Ghahremani M, Nouri K, Whelan JT, Duong J, Tucker S, Diallo JS, Bell JC, Ilkow CS. Implications for SARS-CoV-2 Vaccine Design: Fusion of Spike Glycoprotein Transmembrane Domain to Receptor-Binding Domain Induces Trimerization. *Membranes (Basel)*. 2020 Aug 30;10(9):E215. View article

Rosor AM, Sleigh JN, Groves M, Muntoni F, Reilly MM, Hoogenraad CC, Schiavo G. Loss of BICD2 in muscle drives motor neuron loss in a developmental form of spinal muscular atrophy. *Acta Neuropathol Commun*. 2020 Mar 17;8(1):34. View article

Lester S, Harcourt J, Whitt M, Al-Abdely HM, Midgley CM, Alkhamis AM, Aziz Jokhdar HA, Assiri AM, Tamin A, Thornburg N. Middle East respiratory coronavirus (MERS-CoV) spike (S) protein vesicular stomatitis virus pseudoparticle neutralization assays offer a reliable alternative to the conventional neutralization assay in human seroepidemiological studies. *Access Microbiol*. 2019 Sep 11;1(9):e000057. View article

Yamaguchi H, Honda S, Torii S, Shimizu K, Katoh K, Miyake K, Miyake N, Fujikake N, Sakurai HT, Arakawa S, Shimizu S. Wipi3 is essential for alternative autophagy and its loss causes neurodegeneration. *Nat Commun*. 2020 Oct 20;11(1):5311. View article

Chen DY, Khan N, Close BJ, Goel RK, Blum B, Tavares AH, Kenney D, Conway HL, Ewoldt JK, Kapell S, Chitalia VC, Crossland NA, Chen CS, Kotton DN, Baker SC, Connor JH, Douam F, Emili A, Saeed M. SARS-CoV-2 desensitizes host cells to interferon through inhibition of the JAK-STAT pathway. *bioRxiv [Preprint]*. 2020 Oct 28:2020.10.27.358259. View article

Wei J, Alfajaro MM, DeWeirdt PC, Hanna RE, Lu-Culligan WJ, Cai WL, Strine MS, Zhang SM, Graziano VR, Schmitz CO, Chen JS,

Mankowski MC, Filler RB, Ravindra NG, Gasque V, de Miguel FJ, Patil A, Chen H, Oguntuyo KY, Abriola L, Surovtseva YV, Orchard RC, Lee B, Lindenbach BD, Politi K, van Dijk D, Kadoch C, Simon MD, Yan Q, Doench JG, Wilen CB. Genome-wide CRISPR Screens Reveal Host Factors Critical for SARS-CoV-2 Infection. Cell. 2020 Oct 20;S0092-8674(20)31392-1. View articleVogel AB, et al. BNT162b vaccines protect rhesus macaques from SARS-CoV-2. Nature. 2021 Apr;592(7853):283-289. View articleIf you publish research with this product, please let us know so we can cite your paper.

#### 主要内容

该单克隆抗体与印第安纳血清型的VSV-G蛋白反应。实验：表现出适用于Western印迹，流式细胞术和免疫荧光应用特征的对VSV-G伪型慢病毒（参见以下参考文献）的反应性，可从我们的姐妹公司获得绝对抗体：制造使用具有来自杂交瘤8G5F11VERSULAT is病毒（VSV）的可变区（即特异性）的绝对抗体的重组平台是研究良好的，包膜，阴茎RNA病毒。对于5个蛋白质的VSV基因组编码：N，P，M，G和L.g蛋白（糖蛋白）位于病毒米氏菌表面，并负责病毒附着和渗透。另外，许多慢性载体与来自印第安纳血清型的VSV-G具有假型。从乌沟S. Lyles, Phd, Wake Forest医学院的实验室。

#### 厂牌介绍

### 关于Kerafast Inc.

Kerafast 是一家位于波士顿的试剂公司，其主要使命是为QuanQiu科学界提供易于使用的独特实验室研究工具。我们的产品组合包括细胞系、抗体、小分子、染料等，其中许多在其他地方无法获得。自2011年成立以来，来自[全球190多个机构](#)的研究人员通过我们的在线平台提供了他们的创新试剂，无需通过传统的材料转让协议流程即可快速获取材料。

我们处理提供实验室的所有销售和运输物流，并从每次销售中返还丰厚的特许权使用费。因此，我们帮助提供实验室节省时间和资源，同时为进一步研究提供额外资金。采购科学家可以更轻松地发现和获取其他地方通常无法获得的独特试剂，同时还可以资助其他研究人员的工作。这创建了一个QuanQiu科学家社区，他们贡献和获取Reagent for the Greater Good，以加速他们自己的研究以及整体科学进步。

2018年，Kerafast与[Absolute Antibody](#)合并，后者是一家总部位于英国的公司，其愿景是为所有研究人员提供重组抗体技术。[此次合并](#)将两家公司聚集在一起，共同致力于改善科学界可用的研究工具的选择。

#### 品牌标识



#### 产品关键词

- [kerafast抗体ED2003](#)
- [kerafast 抗体](#)
- [kerafast国内代理商](#)
- [kerafast](#)
- [kerafast 代理](#)
- [kerafast细胞](#)
- [kerafast抗体代理](#)
- [kerafast代理商](#)
- [kerafast品牌](#)
- [kerafast代理](#)
- [kerafast细胞代购kerafast品牌代理](#)
- [kerafast华北代理](#)
- [美国kerafast公司](#)
- [kerafast公司](#)
- [进口kerafast代理](#)
- [kerafast丁香通](#)
- [kerafast中国代理](#)
- [kerafast官网](#)
- [kerafast抗体](#)
- [kerafast专业代理](#)

一键获取大包装优惠报价

#### 选择您的报价场景

- 【我们直接使用】需要优惠报价、大包装规格、货期 -- ---->[报价默认含增值税13%发票；尽量提供货号、规格、需求数量]
- 【需要技术文档】产品说明书、COA、MSDS、手册 -- ---->[默认提供说明书或者COA，特别技术指标要求请下面填入详细描述]
- 【我帮客户找货】需要优惠报价、大包装规格、货期 -- ---->[报价默认含增值税13%发票]
- 【推荐替代产品】需要优惠报价、大包装规格、货期 -- ---->[提供替代产品的价格，默认含增值税13%发票]
- 【我能原厂直采】请只提供代理进口清关服务的报价 -- ---->[适合只需要进口许可证代办服务、清关服务的专业级买家，独立服务]

○ 【其它报价场景】

请输入您的情况与报价要求

报价要求详细描述

贵单位贵姓

接受报价的E-mail

手机扫描二维码阅读本页



可能感兴趣的内容

[NR-19475 土拉弗朗西斯菌亚种.土拉菌,菌株SCHUS4, Gateway?克隆集,在大肠杆菌中重组,板18\(克隆\)](#)

2022-04-01

[NR-3639 甲型流感病毒, Kilbourne F163:](#)

[A/Dunedin/6/1983\(HA,NA\)xA/Puerto Rico/8/1934\(H1N1\), Reassortant X-81\(病毒\)](#)

2022-03-31

[NR-13634 结核分枝杆菌、CDC1551、转座子突变体1106\(MT0951、Rv0924c\)\(突变细菌\)](#)

2022-04-01

[NR-46411 金黄色葡萄球菌, HIP11983\(细菌\)](#)

2022-04-01

[核提取分离培养基, 100ml\(100次反应\)](#)

2021-12-21

[人工模拟合成汗DN 53160-2/BS EN 1811 \(BZ125\) 500ml](#)

2021-12-13

[质谱法, 抗氧化蛋白缓冲液, 100UL](#)

2021-12-21

[NR-781 牛痘病毒, 单克隆抗牛痘\(WR\)A27L, 残基1至110\(腹水\)\(类似于VMC-50\)\(单克隆抗体\)](#)

2022-03-31

[NR-46954 金黄色葡萄球菌亚种.金黄色葡萄球菌、JE2、转座子突变体NE411\(SAUSA300\\_2578\)\(突变细菌\)](#)

2022-04-01

[抗D\[Rho\]抗体, 人73/517](#)

2024-05-19

[NR-41926 肺炎克雷伯菌, BIDMC10\(细菌\)](#)

2022-04-01

[人肾素\[国际标准168/356](#)

2024-05-19

[道康宁PDMS184光学胶灌封胶聚二甲基硅氧烷 1.1KG](#)

2021-12-02

[NR-43499\\_霍氏博德特氏菌, 41130\(细菌\)](#)

2022-04-01

[SARS-CoV-2XF重组分离株101058](#)

2024-05-19

[NR-47932 金黄色葡萄球菌亚种.金黄色葡萄球菌, USA300JE2, 转座子突变体SAUSA300\\_0709\(NE1390\)\(突变细菌\)](#)

2022-04-01

[MRA-105 冈比亚按蚊, M2, 鸡蛋\(矢量\)](#)

2022-04-01

[牛疱疹病毒1型 \(BHV-1 / IBR\) MAb gB-gl IgG2b同种型](#)

2019-05-08

[鼠MET\(V1110i\)表达NIH3T3细胞系, 1个小瓶](#)

2021-12-21

[2024\\_06\\_05\\_Mag](#)

2024-06-03